Celltac chemi CRP 4N

CR-420W, CR-421W

English

General

Intended Purpose

For in vitro diagnostic use only.

Celltac chemi CRP 4N is for in vitro diagnostic use in the automated quantitative measurement of C-reactive protein (CRP) in EDTA treated human whole blood, plasma and serum, as an aid to diagnosis in the detection and evaluation of infection, tissue injury inflammatory disorders and associated diseases. For use with CHM-4100 clinical chemistry analyzer (only when software Ver.03-01 or later is installed) and MEK-1303 automated hematology and clinical chemistry analyzer. Increases in C reactive protein (CRP) values are non-specific for many disease processes and should not be interpreted without a complete clinical evaluation. The Celltac chemi CRP 4N is intended for laboratory professional use only. Testing population are general population patients undergoing evaluation for the associated conditions.

Read the analyzer operator's manual together with this manual before and during use.

- NOTE**
- Nihon Kohden representative will replace a product which proves to be defective until the expiration date shown on the package, provided that the product is used as prescribed by the operating instructions in the operator's manuals.
 - Overall judgement must be performed by the physician, referring to the analysis result, clinical findings and other examination results.
 - Use a test cartridge for a single test. Do not reuse it.
 - Use the product with the specified analyzers only.
 - Read the SDS (Safety Data Sheet) carefully before use. The SDS is available from your Nihon Kohden representative.

Devices Intended for Use in Combination with

- CHM-4100 clinical chemistry analyzer (only when software Ver. 03-01 or later is installed)
- MEK-1303 automated hematology and clinical chemistry analyzer

Materials Provided and Materials Required

Materials Provided

Celltac chemi CRP 4N contains the following reagents used for CRP measurement in a single cartridge.

- CRP diluent
- CRP latex reagent
- CRP hemolysis reagent

Materials Required (Not Provided with the Reagent)

- CR-CAL CRP calibrator
- CRP control (commercially available control materials)

Summary and Explanation

C-reactive protein (CRP) is a type I acute phase reactant, produced by the liver and released into the blood within a few hours after tissue injury, the start of an infection, or other cause of inflammation ¹.

A high level of CRP in the blood is a sign that there may be an inflammatory process occurring in the body. Inflammation itself is not typically a problem, but it can indicate a host of other health concerns, including infection, arthritis, kidney failure, and pancreatitis. CRP can provide prognostic information about risk of future coronary events in apparently healthy subjects.

The reference ranges of serum CRP levels in the Western and Asian populations are widely different Asian population has much lower serum CRP levels (less than one-tenth) than Western populations. ²

¹ Husain, T.M. and Kim, D.H. (2002) C-reactive protein and erythrocyte sedimentation rate in orthopedics. The University of Pennsylvania Orthopaedic Journal, 15, 13-16.

² Kao et al. Serum C-Reactive Protein as a Marker for Wellness Assessment. Annals of Clinical & Laboratory Science, vol 36, no. 2, 2006, 163-169

Analyte or Marker

The Celltac Chemi CRP 4N device is designed to detect C-reactive protein (CRP), a protein made by the liver. CRP levels in the blood increase when there is a condition causing inflammation somewhere in the body. CRP is an acute phase protein, non-specific indicator of inflammation.

Target Treated Population

Testing population are the patients undergoing evaluation for the associated conditions. General population (adults and pediatrics).

Specimen Collection and Preparation

- Use a sample of whole blood, serum or plasma.
- When collecting whole blood or plasma specimen, use the K₂-EDTA, K₃-EDTA or Na₂-EDTA blood collection tube.
- Use the sample within the following period after drawing blood:
 - When storing the sample at 4°C (39°F) or less: 5 days
 - When storing the sample at room temperature (30°C (86°F) or less): 24 hours
- For specimen collection and transportation, follow the instructions for use of the EDTA-blood collection tube or serum-blood collection tube and general recommendations for blood sampling.

Intended Users

For laboratory professional use only, in laboratories with suitable equipment for haematological testing. Qualified personnel, e.g. laboratory technicians trained in haematology analysis techniques, will be able to use according to the operator's manual.

Symbols

The following symbols are used with the measurement kit. The descriptions of each symbol are given in the table below.

Symbol	Description
	Use by
	Lot number
	Catalogue number
	Do not reuse
	Fragile
	Keep away from sunlight
	Temperature limits
	Contains sufficient for <n> tests
	Caution
	Operator's manual; operating instructions
	Model number

Symbol	Description
	Unique device identifier
	In vitro diagnostic medical device
	Exclamation mark ¹
	Manufacturer
	European representative
	The CE mark is a protected conformity mark of the European Union. The four digits after the CE mark indicate the identification number of the Notified Body involved in assessing the product's conformity as a medical device.

Safety Information

-  **WARNING** A warning alerts the user to the possible injury or death associated with the use or misuse of the instrument.
-  **CAUTION** A caution alerts the user to possible injury or problems with the instrument associated with its use or misuse such as instrument malfunction, instrument failure, damage to the instrument, or damage to other property.

Pay attention to all safety information in this operator's manual.

WARNING

Always wear rubber gloves to protect yourself from infection when handling and measuring blood samples and test cartridge. The waste fluid may leak from the test cartridge after use.

CAUTION

- Do not swallow the reagent. If swallowed, rinse the mouth immediately. Do not force vomiting. See a physician.
- If the reagent contacts the eyes, mouth or skin, wash immediately with plenty of water and see a physician.

- NOTE**
- Do not use the product if it is past the expiration date on the package or if the product is stored under unspecified conditions.
 - Store the product inside a polystyrene container and put it in a pharmaceutical refrigerator (2 to 8°C, 36 to 46°F) which allows temperature adjustment. Do not store the product in a domestic refrigerator. This may freeze the reagent.
 - Do not freeze the product. Do not use it if it is frozen.
 - The usage environment of the product is between 15 and 30°C (59 and 86°F).

Hazards Identification

Signal Word

Warning

Hazard Statement

- H317 May cause an allergic skin reaction.

Precautionary Statement Prevention

- P261 Avoid breathing dust/fume/gas/mist/vapours/spray.

¹ The label is affixed in accordance with GHS requirements. For details, refer to "Hazards Identification" (p. 2).

- P272 Contaminated work clothing should not be allowed out of the workplace.
- P280 Wear protective gloves/protective clothing/eye protection/face protection.

Precautionary Statement Response

- P302 + P352 IF ON SKIN: Wash with plenty of water.
- P333 + P313 If skin irritation or rash occurs: Get medical advice/attention.
- P362 + P364 Take off contaminated clothing and wash it before reuse.

Precautionary Statement Disposal

- P501 Dispose of contents/container in accordance with local and national regulations.

2-methylisothiazol-3(2H)-one $\geq 0.0015\%$

Measurement

NOTE: Refer to the analyzer operator's manual for details on measurement.

Measurement Principle

Celltac chemi CRP 4N (CR-420W/421W) uses the Latex agglutination immunoturbidimetric assay method to measure CRP (C-reactive protein) in the blood.


CRP diluent and CRP hemolyzing reagent are added to a whole blood, serum or plasma sample. After hemolysis, CRP latex reagent is added. The CRP in the sample reacts with CRP antibody in the CRP latex reagent, and it agglutinates.

The change in optical absorbance of the agglutinated clumps at a fixed time is measured and this provides the concentration of CRP in the sample.

Procedure

1. Prepare the analyzer. Refer to the analyzer operator's manual.
2. Take a test cartridge out of the refrigerator. Wipe off condensation with a dry cloth.
3. Set the test cartridge into the analyzer. Gently invert the test cartridge several times.
If the agitation is not enough, the result may be inaccurate.

Invert the cartridge more than 3 times.

The analyzer starts reading the QR code on the cartridge, which includes the calibration curve information of the reagent. If an error message appears, refer to the analyzer operator's manual.


4. Measure a sample of whole blood, serum or plasma. When using the whole blood sample, gently invert and mix the sample before use because the red cell components settle easily and become unhomogenous.

The measurement results are shown in mg/dL or mg/L (1 mg/dL = 10 mg/L, 1 mg/L = 0.1 mg/dL). The unit corresponds to the CRP concentration in serum or plasma. Refer to the analyzer operator's manual for details.

- NOTE
- No dirt or damage on the test cartridge
 - No leakage from the test cartridge
 - Before setting the test cartridge, check that the reagent is at the bottom of the wells and that there are no bubbles at the bottom of the wells. If the reagent is at the top of the wells or there are bubbles at the bottom of the wells, gently shake the cartridge to make the reagent go down. Otherwise the test result may be inaccurate.

- Do not use a test cartridge which fell on the floor because the light path surface may be dirty.
- Do not scratch the test cartridge.
- Do not touch the QR code and optical path side of the cartridge.

- In rare cases, depending on the measured sample, the following may occur and cause incorrect measurement.
 - The reagent reacts to an unwanted component in the sample.
 - Interference reaction occurs during measurement.
 - If the measurement value is too high or too low for the patient condition or the measurement result is suspicious, remeasure the sample or measure the sample by another measurement method.

Calibration

Perform calibration in the following cases.

- Calibration is required as a result of regular quality control.
- Calibration is required as a result of quality control performed after maintenance.
- Other cases where calibration is required.

Reagents Used

Use a CRP calibrator in combination with a CR-420W or CR-421W Celltac chemi CRP 4N.

Frequency of Performing Calibration

Determine according to inspections performed by each laboratory.

Procedure

Refer to the analyzer operator’s manual or service manual.

Quality Control

Measurement of a control is recommended for statistical quality control of the measurement accuracy of the analyzer. Also, it is recommended to use the average and upper and lower limits set by each laboratory to measure the control.

Frequency of Performing Quality Control

Determine according to the procedures established by each laboratory.

Control Used

Use commercially available control materials.

Procedure

Refer to the analyzer operator’s manual.

Variable Factors

Variable factors that should be considered when setting appropriate upper and lower limits for quality control include the following.

- Calibrator lots
- Control lots etc.
- Reagent lots
- Analyzer used for the control measurements
- Measurement frequency

Technical Information

Metrological Traceability of Values Assigned to the Calibrator

For CRP, the internationally recognized reference method is: IRMM: Institute for Reference Materials and Measurements, European Commission Joint Research Centre Certificate of Analysis ERMDA470/ IFCC and ERM-DA474/IFCC.

Quality Control, Analytical Performance Characteristics

Sensitivity Test

Variation in optical absorbance is between 0.05 and 0.5 when a CRP calibrator of 4.0 mg/dL is used as the sample.

Accuracy Test

Test result is ±15% of the known concentration when a control sample of a known concentration is used as the sample.

Repeatability Test

CV (coefficient of variation) of test result is within 10% when a control sample of a known concentration is assayed 5 times.

Measurement Range (Measuring Interval)

The lower limit of the measuring interval is defined by the limit of quantitation (LoQ). The upper limit of the measuring interval was established through a linearity study.

- Serum or plasma (hematocrit 0%): 1 to 200 mg/L ¹
- Whole blood (hematocrit 40%): 1 to 330 mg/L ¹

¹ Confirmed on the MEK-1303 and CHM-4100

Expected Values

Expected values in normal populations/Reference range in serum and plasma

≤ 5 mg/L

Expected values in affected populations

Diagnostic sensitivity, diagnostic specificity, positive predictive value and negative predictive value are not commonly reported for CRP assays. This is largely attributed to the fact that CRP is used in conjunction with other indicators for the intended clinical applications and patient treatment decision-making. Similarly, traditional likelihood ratio is not commonly reported for CRP.

Literature: Guidance for Industry and FDA Staff: Review Criteria for Assessment of C-Reactive Protein (CRP), High Sensitivity C-Reactive Protein (hsCRP) and Cardiac C-Reactive Protein (cCRP) Assays

Batch Variation

Three levels of control was five repeated in the three days. Between lot %CV was below.

Model	Control Level 1	Control Level 2	Control Level 3
MEK-1303	1.6%	1.9%	2.1%
CHM-4100	0.0%	0.4%	0.1%

Composition

- CRP diluent: Buffer solution
- CRP latex reagent: Latex particles coated with rabbit anti-human CRP polyclonal antibodies: 0.32 ±0.08 mg (per 1 test cartridge), preservative agent: Proclin 950
- CRP hemolysis reagent: Buffer solution containing surfactant, preservative agent: Proclin 950

Sterilization Method

The CR-420W or CR-421W Celltac chemi CRP 4N is not intended to be sterilized or kept in a sterile environment.

Interfering Substances or Limitations

The following interfering substances have been confirmed to have no effect on test results below the listed concentrations.

Substance	Max. Concentration
Rheumatoid factor	774 IU/mL
Bilirubin C (conjugated type)	30.8 mg/dL
Bilirubin F (free type)	27.6 mg/dL
Chyle	2,344 FTU
Ascorbic acid	322 mg/dL

Analytical Performance Characteristics

Values shown in this section are representative performance data of the Celltac Chemi CRP 4N and may differ from values acquired in individual facilities.

Analytical Sensitivity

To determine Limit of Detection (LoD), Limit of Blank (LoB), and Limit of Quantification (LoQ) of Celltac chemi CRP 4N on MEK-1303 and CHM-4100.

Reference guidance: CLSI EP17-A2

Model	LoB	LoD	LoQ
MEK-1303	0.003 mg/dL	0.012 mg/dL	0.051 mg/dL
CHM-4100	0.016 mg/dL	0.027 mg/dL	0.073 mg/dL

Precision

One lot of Celltac Chemi CRP 4N was used to evaluate precision. Three levels of Control was repeated twice in the morning and evening for 25 days.

Liquichek Elevated CRP Control: Bio-Rad Laboratories Inc.

Reference guidance: CLSI EP05-A3, CLSI H26-A2

Model	Sample		Mean	n	Repeatability		Within-laboratory Precision	
					SD	%CV	SD	%CV
MEK-1303	Hematology Control + CRP	MK-3CL	0.504 mg/dL	100	0.017	3.4%	0.030	6.0%
		MK-3CN	2.639 mg/dL	100	0.053	2.0%	0.095	3.6%
		MK-3CH	6.150 mg/dL	100	0.147	2.4%	0.197	3.2%
CHM-4100	Liquichek Elevated CRP Control	Level 1	1.239 mg/dL	100	0.0235	1.9%	0.0286	2.3%
		Level 2	8.023 mg/dL	100	0.1132	1.4%	0.1532	1.9%
		Level 3	16.729 mg/dL	100	0.2373	1.4%	0.3265	2.0%

Accuracy, Trueness

For Celltac chemi CRP 4N assay trueness is demonstrated through traceability to the certified reference material and comparison study of the device to a commercially available similar device presented.

Actual measurements, their mean, and 95%CI shall be within the certified value ± uncertainty.

The IFCC blood serum CRP international standard (ERM-DA474/IFCC) is measured 5 times in the serum mode, using 3 different batches of CR-420W assay.

Model	ERM-DA474 Certified Value (mg/dL)	Reagent Lot		
		1	2	3
MEK-1303	4.12 +/-0.25	4.13-4.21 mg/dL	4.20-4.34 mg/dL	4.13-4.19 mg/dL

The IFCC blood serum CRP international standard (ERM-DA474/IFCC) is measured 5 times in the serum mode, using 3 combinations of reagent and instrument (CHM-4100).

Model	ERM-DA474 Certified Value (mg/dL)	Combination		
		1	2	3
CHM-4100	4.12 +/-0.25	4.15-4.21 mg/dL	4.15-4.21 mg/dL	4.06-4.08 mg/dL

Linearity

Samples were measured in quadruplicate using the MEK-1303 and CHM-4100 analyzers to confirm linearity across the measuring range.

Reference guidance: CLSI EP06-A

This kit measuring interval obtained was below:

Model	Linearity Measuring Interval
MEK-1303	0.003-20.410 mg/dL
CHM-4100	0.008-20.233 mg/dL

Correlation

Correlation coefficient and slope calculated by least squares method

Reference guidance: CLSI EP09-A3

Model	n	Correlation Coefficient	Intercept	Slope
MEK-1303 v.s. Other company	131	0.998	0.006	0.9656
CHM-4100 v.s. Other company	131	0.999	-0.011	0.9246

High Dose Hook Effect

Samples with a CRP concentration of about 30 mg/dL and about 40 mg/dL were measured three times using the CR-420W/421W on the MEK-1303.

Hook effect is not confirmed until CRP concentrations of 40 mg/dL

- CLSI. *Evaluation Of Detection Capability For Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. CLSI document EP17-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2012.
- CLSI. *Validation, Verification, and Quality Assurance of Automated Hematology Analyzers; Approved Standard—Second Edition*. CLSI document H26-A2. Wayne, PA: Clinical and Laboratory
- CLSI. *Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline—Third Edition*. CLSI document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
- CLSI. *Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline*. CLSI document EP06-A. Wayne, PA: Clinical and Laboratory Standards Institute; 2003.
- CLSI. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. CLSI document EP09-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.

Mathematical Approach upon which the Calculation of the Analytical Result is Made

Calibration information is determined at the factory for each reagent lot. The calibration information is recorded in the QR code, and the analyzer reads the information before measurement and calculates the CRP value.

Environmental Conditions

Storage and Transport Environment

Temperature: 2 to 8°C (36 to 46°F)

Usage Environment

Temperature: 15 to 30°C (59 to 86°F)


Expiration Date

The expiration date is shown on the package.


Lifetime: 18 months

Package and Catalog Number


CRP diluent, CRP latex reagent and CRP hemolysis reagent are contained in a single test cartridge.

CR-420W

CR-421W

Test cartridge

Model	Qty	Catalog Number
CR-420W	50 (25 test cartridges × two boxes)	CR-420W
CR-421W	25	CR-421W

Disposal

⚠ WARNING

Dispose of the test cartridge according to your local laws and your facility's guidelines (including incineration, melt treatment, sterilization, disinfection and request for waste disposal) for **disposing of infectious medical waste**. Otherwise, it may affect the environment. If there is a possibility that the product may have been contaminated with infection, it may cause infection.

When disposing of the test cartridges, such as when the expiration date is past, follow the instructions on the SDS of the Celltac chemi CRP 4N.

Revision History

Edition	Date	Details	Code Number
1st Edition	24 Jan 2018	Initial issue	0654-905395
5th Edition	12 Dec 2022	Symbols added	0614-907956B
6th Edition	08 Jan 2025	Corrections	0614-907956C

- NOTE
- The code number of this manual was changed from 0654-905395B to 0614-907956A when the manual was updated from 3rd Edition to 4th Edition.
 - Changes made in the most recent edition are indicated by a bar in the left margin of each page.

Trademark

“QR Code” is registered trademark of DENSO WAVE INCORPORATED.

Note for users in the territory of the EEA and Switzerland:
Any serious incident that has occurred in relation to the device must be reported to the European Representative designated by the manufacturer and the Competent Authority of the Member State of the EEA and Switzerland in which the user and/or patient is established.

Copyright Notice

The entire contents of this manual are copyrighted by Nihon Kohden. All rights are reserved.

Manufacturer

NIHON KOHDEN CORPORATION

1-31-4 Nishi-chiai, Shinjuku-ku,
Tokyo 161-8560, Japan
Phone +81 3-5996-8041
<https://www.nihonkohden.com/>

European Representative

NIHON KOHDEN EUROPE GmbH

Raiffeisenstrasse 10, 61191 Rosbach, Germany
Phone +49 6003-827-0 Fax +49 6003-827-599

UK Responsible Person

NIHON KOHDEN UK LTD.

Unit 3, Heyworth Business Park,
Old Portsmouth Road, Peasmarsh,
Guildford, Surrey, GU3 1AF, UK
Phone +44 14-8333-1328

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- CLSI. *Evaluation Of Detection Capability For Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. CLSI-Dokument EP17-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2012.
- CLSI. *Validation, Verification, and Quality Assurance of Automated Hematology Analyzers; Approved Standard—Second Edition*. CLSI-Dokument H26-A2. Wayne, PA: Clinical and Laboratory
- CLSI. *Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline—Third Edition*. CLSI-Dokument EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
- CLSI. *Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline*. CLSI document EP06-A. Wayne, PA: Clinical and Laboratory Standards Institute; 2003.
- CLSI. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. CLSI-Dokument EP09-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.

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- CLSI. *Evaluation Of Detection Capability For Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. CLSI document EP17-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2012.
- CLSI. *Validation, Verification, and Quality Assurance of Automated Hematology Analyzers; Approved Standard—Second Edition*. CLSI document H26-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
- CLSI. *Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline—Third Edition*. CLSI document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
- CLSI. *Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline*. CLSI document EP06-A. Wayne, PA: Clinical and Laboratory Standards Institute; 2003.
- CLSI. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. CLSI document EP09-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.

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- CLSI. *Evaluation Of Detection Capability For Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. CLSI document EP17-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2012.
- CLSI. *Validation, Verification, and Quality Assurance of Automated Hematology Analyzers; Approved Standard—Second Edition*. CLSI document H26-A2. Wayne, PA: Clinical and Laboratory
- CLSI. *Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline—Third Edition*. CLSI document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
- CLSI. *Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline*. CLSI document EP06-A. Wayne, PA: Clinical and Laboratory Standards Institute; 2003.

- CLSI. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. CLSI document EP09-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.

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- CLSI. *Evaluation Of Detection Capability For Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. CLSI document EP17-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2012.
- CLSI. *Validation, Verification, and Quality Assurance of Automated Hematology Analyzers; Approved Standard—Second Edition*. CLSI document H26-A2. Wayne, PA: Clinical and Laboratory
- CLSI. *Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline—Third Edition*. CLSI document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.

- CLSI. *Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline*. CLSI document EP06-A. Wayne, PA: Clinical and Laboratory Standards Institute; 2003.
- CLSI. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. CLSI document EP09-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.

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Efekat kuke nije potvrđen sve do koncentracija CRP-a od 40 mg/dl

- CLSI. *Evaluation Of Detection Capability For Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. CLSI document EP17-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2012.
- CLSI. *Validation, Verification, and Quality Assurance of Automated Hematology Analyzers; Approved Standard—Second Edition*. CLSI document H26-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2012.
- CLSI. *Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline—Third Edition*. CLSI document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
- CLSI. *Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline*. CLSI document EP06-A. Wayne, PA: Clinical and Laboratory Standards Institute; 2003.
- CLSI. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. CLSI document EP09-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.

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- CLSI. *Evaluation Of Detection Capability For Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. Documento EP17-A2 del CLSI. Wayne, PA: Clinical and Laboratory Standards Institute; 2012.
- CLSI. *Validation, Verification, and Quality Assurance of Automated Hematology Analyzers; Approved Standard—Second Edition*. Documento H26-A2 del CLSI. Wayne, PA: Clinical and Laboratory
- CLSI. *Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline—Third Edition*. Documento EP05-A3 del CLSI. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.

- CLSI. *Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline*. Documento EP06-A del CLSI. Wayne, PA: Clinical and Laboratory Standards Institute; 2003.
- CLSI. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. Documento EP09-A3 del CLSI. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.

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