CRP Calibrator
CR-CAL

English

General

Intended Purpose

For in vitro diagnostic use only.

CRP calibrator is for in vitro diagnostic use in the calibrations of the CR-420W and CR-421W CRP assay.

Devices Intended for Use in Combination with

- CHM-4100 clinical chemistry analyzer (only when software Ver. 03-01 or later is installed)
- MEK-1303 automated hematology and clinical chemistry analyzer

NOTE • The specifications and equipment requirements can be found in the operator's manual for the above mentioned Nihon Kohden in vitro diagnostic devices.

- Read the SDS (Safety Data Sheet) carefully before use. The SDS is available from your Nihon Kohden representative.

Materials Provided and Materials Required

Materials Provided

CR-CAL CRP calibrator

Materials Required (Not Provided with the Calibrator)

CR-420W or CR-421W Celltac chemi CRP 4N

Intended Users

For laboratory professional use only, in laboratories with suitable equipment for haematological testing. Qualified personnel, e.g. laboratory technicians trained in haematology

analysis techniques, will be able to use according to the operator's manual.

Symbols

The following symbols are used with this calibrator. The descriptions of each symbol are given in the table below.

Symbols	Description	Symbols	Description
	Caution		Biological risks
	Use by		Operator's manual; operating instructions
	Lot number		Keep away from sunlight
	Catalogue number		Temperature limits
	Manufacturer		This way up
	European representative		Fragile
	The CE mark is a protected conformity mark of the European Union. The four digits after the CE mark indicate the identification number of the Notified Body involved in assessing the product's conformity as a medical device.		Keep away from rain
	Model number		Stacking limit by number ("n" is the limiting number)
	Unique device identifier		
	In vitro diagnostic medical device		

Safety Information

⚠ WARNING A warning alerts the user to possible injury or death associated with the use or misuse of the instrument.

Pay attention to all safety information in this operator's manual.

⚠ WARNING

Do not swallow the calibrator. It may be infectious.

⚠ WARNING

When handling the calibrator, wear rubber gloves to protect yourself from infection.

⚠ WARNING

POTENTIALLY BIOHAZARDOUS MATERIAL. This calibrator is potential biohazardous material. The CRP calibrator is intended solely for in vitro diagnostic use by trained personnel. The calibrator contains human-sourced and/or potentially infectious components. Components from human donors used in preparation of the calibrator were tested for the presence of the antibodies to human immunodeficiency virus (HIV-1/HIV-2) and Hepatitis C virus (HCV) as well as for Hepatitis B virus surface antigen and found to be negative. Because no test method can offer complete assurance that infectious agents are absent, this material should be handled as potentially infectious. When handling or disposing of vials follow precautions for patient specimens as specified in the OSHA Bloodborne Pathogen Rule (29 CFR Part 1910, 1030) or other equivalent biosafety procedures.

- Do not use the calibrator after the expiration date.
- Store the calibrator in a pharmaceutical refrigerator.
- Store the calibrator at temperatures between 2 and 8°C (36 and 46°F).
- Do not freeze the calibrator vial.

- If the cap of the calibrator becomes dirty, wipe with a non-woven towel such as Kimwipes.

Hazards Identification

Supplemental Hazard Information

EUH032 Contact with acids liberates very toxic gas.

Sodium azide <0.1%

Using the Calibrator

Perform calibration in the following cases.

- Calibration is required as a result of regular quality control.
- Calibration is required as a result of quality control performed after maintenance.
- Other cases where calibration is required.

Reagents Used

Use the CRP calibrator in combination with a CR-420W or CR-421W Celltac chemi CRP 4N.

Frequency of Performing Calibration

Determine according to inspections performed by each laboratory.

Measurement Principle


Refer to the analyzer operator's manual.

Procedure

NOTE • Turbidity, discoloration of the calibrator or unacceptable results may indicate deterioration. Do not use the calibrator if deterioration is suspected.

- Confirm that the lot number of the calibrator matches the lot number on the calibration data label.


1. Take the calibrator out of the medical refrigerator and allow it to warm to room temperature for 15 minutes before mixing.

2. Enter “Assay value” to the analyzer, refer to calibration data label.
3. Invert and mix the calibrator. (more than 8 times)

NOTE

- Do not mix the calibrator using a mixer.
- Immediately before measuring the calibrator, gently turn it upside down.
- Incomplete mixing of the calibrator affects the measurement result.

4. Open the cap and aspirate the calibrator into the analyzer.

NOTE: Insert the sampling nozzle all the way to the bottom of the vial so that the tip is below the liquid surface.

5. Close the cap.

NOTE

- If aspirating the calibrator into the analyzer again, invert and mix with the cap closed.
- Use the calibrator within the expiration date and dispose of it after the expiration date for disposal of infectious medical waste.
- If there are unused vial, restore the calibrator at temperatures between 2 and 8°C (36 and 46°F)
- Follow the instructions in the analyzer operator's manual, or the section titled “Calibration” in the analyzer service manual.”

Technical Information

Metrological Traceability of Values Assigned to the Calibrator

For CRP, the internationally recognized reference method is:

IRMM: Institute for Reference Materials and Measurements, European Commission Joint Research Centre Certificate of Analysis ERM-DA470/IFCC and ERM-DA474/IFCC.

Composition

The calibrator contains fibrin-removed human plasma, human-sourced components, stabilizers and preservatives.

Expected Values

For calibration values, refer to the assay data printed on the calibration data label.

Sterilization Method

The CRP calibrator is not intended to be sterilized or kept in a sterile environment.

Storage and Transport Environment

Transport and store at temperatures between 2 and 8°C (36 to 46°F)

Expiration Date

Shown on the container label.

Shelf Life Date After Opening the Package

Opened calibrator should be used within the day it has been opened.

Package and Catalog Number

Model	Qty	Catalog Number
CR-CAL	0.25 mL vial × 3	CR-CAL

Disposal

⚠ WARNING

Dispose of the calibrator according to your local laws and your facility's guidelines (including incineration, melt treatment, sterilization, disinfection and request for waste disposal) for disposing of infectious medical waste. Otherwise, it may affect the environment. If there is a possibility that the product may have been contaminated with infection, it may cause infection.

When disposing of the calibrator, such as when the expiration date is past, follow the instructions on the SDS of the calibrator.

Revision History

Edition	Date	Details	Code Number
1st Edition	24 Jan 2018	Initial issue	0604-912879
4th Edition	12 Dec 2022	Symbols added	0614-907958B
5th Edition	08 Jan 2025	Corrections	0614-907958C

- NOTE
- The code number of this manual was changed from 0604-912879A to 0614-907958A when the manual was updated from 2nd Edition to 3rd Edition.
 - Changes made in the most recent edition are indicated by a bar in the left margin of each page.

Note for users in the territory of the EEA and Switzerland:
Any serious incident that has occurred in relation to the device must be reported to the European Representative designated by the manufacturer and the Competent Authority of the Member State of the EEA and Switzerland in which the user and/or patient is established.

 Manufacturer
NIHON KOHDEN CORPORATION
1-31-4 Nishiochiai, Shinjuku-ku,
Tokyo 161-8560, Japan
Phone +81 3-5996-8041
<https://www.nihonkohden.com/>

EC REP European Representative
NIHON KOHDEN EUROPE GmbH
Raiffeisenstrasse 10, 61191 Rosbach, Germany
Phone +49 6003-827-0 Fax +49 6003-827-599
UK Responsible Person
NIHON KOHDEN UK LTD.
Unit 3, Heyworth Business Park,
Old Portsmouth Road, Peasmarsh,
Guildford, Surrey, GU3 1AF, UK
Phone +44 14-8333-1328

Copyright Notice
The entire contents of this manual are copyrighted by Nihon Kohden. All rights are reserved.

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