HA-410V CAL

YZ-005B1

English

General

Intended Purpose

For in vitro diagnostic use only.

HA-410V CAL is for in vitro diagnostic use in the calibration of the HA-420W and HA-421W HbA1c assay. HA-410V CAL is an HbA1c calibrator to be primarily used by domestic and overseas service personnel for calibration of Clinical Chemistry Analyzer CHM-4100 and the Automated Hematology and Clinical Chemistry Analyzer MEK-1303 for regular maintenance and repairs.

Devices Intended for Use in Combination with

- CHM-4100 clinical chemistry analyzer
- MEK-1303 automated hematology and clinical chemistry analyzer

- NOTE**
- The specifications and equipment requirements can be found in the operator's manual for the abovementioned Nihon Kohden in vitro diagnostic devices.
 - Read the SDS (Safety Data Sheet) carefully before use. The SDS is available from your Nihon Kohden representative.
 - The SSP (Summary of Safety and Performance) for the HA-410V CAL can be obtained from the following URL:
<https://ec.europa.eu/tools/eudamed>

Materials Provided and Materials Required

Materials Provided

YZ-005B1 HA-410V CAL (Three vials of lyophilised calibrator are supplied in each pack.)

Materials Required (Not Provided with the Reagent)

- HA-420W or HA-421W Celltac chemi HbA1c N
- Distilled water
- Micro-pipette or syringe that can aspirate 1.5 mL of fluid (locally purchased)

Intended Users

HA-410V CAL is an HbA1c calibrator to be primarily used by domestic and overseas service personnel for calibration of Clinical Chemistry Analyzer CHM-4100 and the Automated Hematology and Clinical Chemistry Analyzer MEK-1303 for regular maintenance and repairs.

Symbols

The following symbols are used with this calibrator. The descriptions of each symbol are given in the table below.

Symbols	Description
	Caution
	Use by
	Lot number
	Catalogue number
	Manufacturer

Symbols	Description
	Operator's manual; operating instructions
	Model number
	Unique device identifier
	In vitro diagnostic medical device
	Biological risks

Symbols	Description
	Vial
	Temperature limits
	Keep away from sunlight
	This way up
	Fragile
	Keep away from rain
	Stacking limit by number ("n" is the limiting number)

Symbols	Description
	European representative
	The CE mark is a protected conformity mark of the European Union. The four digits after the CE mark indicate the identification number of the Notified Body involved in assessing the product's conformity as a medical device.

Safety Information

⚠ WARNING A warning alerts the user to possible injury or death associated with the use or misuse of the instrument.

Pay attention to all safety information in this operator's manual.

⚠ WARNING

Do not swallow the calibrator. It may be infectious.

⚠ WARNING

When handling the calibrator, wear rubber gloves to protect yourself from infection.

⚠ WARNING

POTENTIALLY BIOHAZARDOUS MATERIAL. This calibrator is potential biohazardous material. The calibrator is intended solely for in vitro diagnostic use by trained personnel. The calibrator contains human-sourced and/or potentially infectious components. Components from human donors used in preparation of the calibrator were tested for the presence of the antibodies to human immunodeficiency virus (HIV-1/HIV-2) and Hepatitis C virus (HCV) as well as for Hepatitis B virus surface antigen and found to be negative. Because no test method can offer complete assurance that infectious agents are absent, this material should be handled as potentially infectious. When handling or disposing of vials follow precautions for patient specimens as specified in the OSHA Bloodborne Pathogen Rule (29 CFR Part 1910, 1030) or other equivalent biosafety procedures.

- Do not use the calibrator after the expiration date.
- If microorganisms get into the calibrator or the calibrator vial is turbid, discard the calibrator.
- The quality control program cannot be performed on this calibrator.
- Store the calibrator in a pharmaceutical refrigerator.
- Allow the calibrator to warm to temperatures between 15 and 30°C (59 and 86°F) before use.
- Invert and mix the calibrator before use.
- The assay value may vary for each vial. Check the value after every measurement and enter it on the analyzer.

Using the Calibrator

Perform calibration in the following cases.

- Calibration is required as a result of regular quality control.
- Calibration is required as a result of quality control performed after maintenance.

- Other cases where calibration is required.

Reagents Used

Use the HA-410V CAL in combination with an HA-420W or HA-421W Celltac chemi HbA1c N.

Frequency of Performing Calibration

Determine according to inspections performed by each laboratory.

Measurement Principle


Refer to the analyzer service manual.

Procedure


NOTE • Turbidity, discoloration of the calibrator or unacceptable results may indicate deterioration. Do not use the calibrator if deterioration is suspected.

- Confirm that the lot number of the calibrator matches the lot number on the calibration data label.


1. Take the calibrator out of the refrigerator.
2. Leave the calibrator for about **10 minutes** to warm it to room temperature (15 to 30°C, 59 to 86°F).

3. Enter "Assay value" to the analyzer, refer to calibration data label.
4. Open the lid of the calibrator vial.
5. Aspirate 1.5 mL of distilled water using the micro-pipette or syringe, and dispense the water into the calibrator vial.

6. Close the lid of the calibrator vial and leave the calibrator for 10 minutes.
7. After 10 minutes have passed, gently invert the calibrator vial 20 times so the solution inside the vial is thoroughly mixed. Be careful not to allow bubbles to form while inverting.

NOTE: Insert the sampling nozzle all the way to the bottom of the vial so that the tip is below the liquid surface.

8. Aspirate the calibrator into the analyzer.
9. Close the cap.

NOTE • If aspirating the calibrator into the analyzer again, invert and mix with the cap closed.

- Use the calibrator within the expiration date and dispose of it after the expiration date for disposal of infectious medical waste.
- When performing the calibration five times, be careful not to take too much time and avoid performing calibration under significant temperature change.
- Invert and mix the calibrator before use.
- Write down the calibrator vial opened date on the label on the calibrator and put the calibrator back to the

refrigerator (2 to 8°C (36 to 46°F)) immediately after the measurement.

- Do not mix the calibrator with another calibrator with a different lot or another reagent.
- Do not reuse the calibrator. Do not add the reagent to the calibrator container.
- If there are unused vial, restore the calibrator at temperatures between 2 and 8°C (36 and 46°F)
- You can use the calibrator up to four hours in room temperature. If you need to use the calibrator more than four hours, put the calibrator in a refrigerator or benchtop cooler during measurement. (If it is used in a cold place, the calibrator can be used more than four hours.)
- Use a new calibrator when calibrating more than one chemistry analyzer or when the calibrator runs out before calibration is finished.

Technical Information

Metrological traceability of values assigned to the calibrator

Assigned values of HbA1c traceable to an internationally recognized reference material: JCCRM-411

The measurement results are reported in %NGSP. This unit corresponds to HbA1c concentration as a percentage of total hemoglobin. The following equation shows the relationship between NGSP and IFCC.

$$\text{NGSP value (\%)} = 0.0915 \times \text{IFCC value (mmol/mol)} + 2.15$$

Results are traceable to HbA1c calibrator for NGSP certification (JCCAL)/Certified Reference Material for Measurement of HbA1c (JCCRM).

The following performance specification is identified for the HbA1c calibrator:

$$\text{HbA1c concentration} = 9.1\text{-}9.9\%$$

Analytical Performance Characteristics

Precision: One lot of Celltac Chemi HbA1c

N (HA-420W/421W) was used to evaluate precision. One lot of HA-410V CAL(YZ-005B1) was repeated five times in the morning and evening for 10 days.

Reference guidance: CLSI EP05-A3

Model		MEK-1303	CHM-4100
Sample		HA-410V CAL (YZ-005B1)	
Mean (%)		9.752	9.559
n		100	100
Repeatability	SD	0.065	0.079
	CV(%)	0.7	0.8
Within-Laboratory Precision	SD	0.096	0.110
	CV(%)	1.0	1.2

CLSI. *Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline—Third Edition*. CLSI document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.

Sterilization Method

The HbA1c calibrator is not intended to be sterilized or kept in a sterile environment.

Composition

Red lyophilized reagent (human whole blood based lyophilized reagent [HbA1c concentration (9.1-9.9%)])

Storage and Transport Environment

2 to 8°C (36 to 46°F)

Expiration Date

Shown on the container label.

Shelf Life Date After Opening the Package

Calibrator should be used within the same day it is thawed.

Once reconstituted, use within the day while stored at 2 to 8°C or within 4 hours when stored at room temperature.

Package and Catalog Number

Model	Qty	Catalog Number
YZ-005B1	3 vials	T458

Revision History

Edition	Date	Details	Code Number
1st Edition	09 Aug 2016	Initial issue	0604-912509
6th Edition	12 Dec 2022	Symbols added	0614-907959B
7th Edition	08 Jan 2025	Corrections	0614-907959C

Disposal

 **WARNING**

Dispose of the calibrator according to your local laws and your facility’s guidelines (including incineration, melt treatment, sterilization, disinfection and request for waste disposal) for disposing of infectious medical waste. Otherwise, it may affect the environment. If there is a possibility that the product may have been contaminated with infection, it may cause infection.

When disposing of the calibrator, such as when the expiration date is past, follow the instructions on the SDS of the calibrator.

- NOTE
- The code number of this manual was changed from 0604-912509E to 0614-907959A when the manual was updated from 4th Edition to 5th Edition.
 - Changes made in the most recent edition are indicated by a bar in the left margin of each page.

Note for users in the territory of the EEA and Switzerland:
Any serious incident that has occurred in relation to the device must be reported to the European Representative designated by the manufacturer and the Competent Authority of the Member State of the EEA and Switzerland in which the user and/or patient is established.

 Manufacturer
NIHON KOHDEN CORPORATION
1-31-4 Nishiochiai, Shinjuku-ku,
Tokyo 161-8560, Japan
Phone +81 3-5996-8041
<https://www.nihonkohden.com/>

EC REP European Representative
NIHON KOHDEN EUROPE GmbH
Raiffeisenstrasse 10, 61191 Rosbach, Germany
Phone +49 6003-827-0 Fax +49 6003-827-599

UK Responsible Person
NIHON KOHDEN UK LTD.
Unit 3, Heyworth Business Park,
Old Portsmouth Road, Peasmarsh,
Guildford, Surrey, GU3 1AF, UK
Phone +44 14-8333-1328

Copyright Notice
The entire contents of this manual are copyrighted by Nihon Kohden. All rights are reserved.

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