HEMATOLOGY CONTROL for WBC 5 part differential MEK-5DL, MEK-5DN, MEK-5DH

English

General

Intended Purpose

For in vitro diagnostic use only.

HEMATOLOGY CONTROL for WBC 5 part differential is for in vitro diagnostic use as a value assigned whole blood control designed to monitor the values of Nihon Kohden hematology analyzers which can measure WBC 5 part differential. HEMATOLOGY CONTROL for WBC 5 part differential is not used for diagnosis or assistance in diagnosis, prognosis, prediction of disease or patient's physiological condition. Measuring this HEMATOLOGY CONTROL for WBC 5 part differential ensures that the device is functioning according to its intended use.

Control parameters are below:

- 1) WBC: White Blood Cell Count
- 2) NE%: Neutrophil Percent
- 3) LY%: Lymphocyte Percent
- 4) MO%: Monocyte Percent
- 5) EO%: Eosinophil Percent
- 6) BA%: Basophil Percent
- 7) NE: Neutrophil
- 8) LY: Lymphocyte
- 9) MO: Monocyte
- 10) EO: Eosinophil
- 11) BA: Basophil
- 12) RBC: Red Blood Cell Count
- 13) HGB: Hemoglobin Concentration
- 14) HCT: Hematocrit Percent
- 15) MCV: Mean Corpuscular Volume
- 16) MCH: Mean Corpuscular Hemoglobin
- 17) MCHC: Mean Corpuscular Hemoglobin Concentration
- 18) RDW: Red Blood Cell Distribution Width
- 19) PLT: Platelet Count
- 20) MPV: Mean Platelet Volume
- 21) IG%: Immature Granulocyte Percent ¹
- 22) IG: Immature Granulocyte ¹

¹ MEK-9200 only

Only software versions 02-01 and later can be used

Read the analyzer operator's manual together with this manual before and during use. Also refer to the assay sheet available on the Nihon Kohden website (See the "Materials Provided" section).

- NOTE
- Use the control with the specified analyzers only.
 - Read the SDS (Safety Data Sheet) carefully before use. The SDS is available from your Nihon Kohden representative and Nihon Kohden website (<https://www.nihonkohden.com>).

Devices Intended for Use in Combination with

MEK-7222, MEK-7300, MEK-8222, MEK-9100, MEK-9200

Materials Provided and Materials Required


Materials Provided

- MEK-5DL, MEK-5DN, MEK-5DH
HEMATOLOGY CONTROL for WBC 5 part differential
- Assay sheet ¹

¹ How to obtain assay sheets

NOTE: Assay sheets are not included. Follow the steps below to obtain assay sheets.

1. Visit our website below to download assay sheets for the control.
<https://www.nihonkohden.com/technology/assaysheet/assaysheet.html>
2. When printing, please use a laser printer.
Set the print quality to 1200×1200 dpi or higher on the printer.

Materials Required (Not Provided with the Control)

- Protective gloves
- Protective glasses
- Waste container for biohazardous materials

Intended Users

For laboratory professional use only, in laboratories with suitable equipment for hematological testing. Qualified personnel, e.g. laboratory technicians trained in hematology analysis techniques, will be able to use according to the operator's manual.

Symbols

The following symbols are used with the control. The descriptions of each symbol are given in the table below.

Symbol	Description	Symbol	Description
	Caution		Biological risks
	Operator's manual; operating instructions		This way up
	Catalogue number		Control/ Concentration level
	Lot number		Open tube stability
	Model number		Unique Device Identifier
	Use by		The CE mark is a protected conformity mark of the European Union. The four digits after the CE mark indicate the identification number of the Notified Body involved in assessing the product's conformity as a medical device.
	In vitro diagnostic medical device		
	Manufacturer		
	Authorized representative in the European Community/ European Union		
	Temperature limits		

Safety Information

WARNING A warning alerts the user to the possible injury or death associated with the use or misuse of the instrument.

Pay attention to all safety information in this operator's manual.

WARNING

- POTENTIALLY BIOHAZARDOUS MATERIAL.** For in vitro diagnostic use. Each human donor/unit used in the preparation of the control has been tested by an FDA licensed method/test and found to be negative or non-reactive for the presence of HBsAg, Anti-HCV, NAT testing for HIV-1, HCV (RNA) and HIV-1/2. Each unit is also negative by a serological test for Syphilis (RPR or STS). Because no test method can offer complete assurance that infectious agents are absent, this material should be handled as potentially infectious. When handling or disposing of tubes follow precautions for patient specimens as specified in the OSHA Bloodborne Pathogen Rule (29 CFR Part 1910, 1030) or other equivalent biosafety procedures.
- Wear protective equipment such as disposable gloves when handling the control.
- Do not swallow the control.

NOTE • Do not use the control if it is past the expiration date on the control or if the control is stored under unspecified conditions.

- Store the control in a pharmaceutical refrigerator.
- Store the control at temperatures between 2 and 8°C (36 and 46°F).
- Do not freeze the control.
- Accurate values cannot be obtained if the control is used in a non-mechanical hematological testing.

Using the Control

NOTE • Refer to the analyzer operator's manual for details on measurement.

- Do not use the control if deterioration is suspected. After mixing, the control should be similar in appearance to fresh whole blood. In unmixed tubes, the supernatant may appear cloudy and reddish; this is normal and does not indicate deterioration. Other discoloration, very dark red supernatant or unacceptable results may indicate deterioration.
- If the control is stored out of the specified temperature range, LY% and LY values may become high.
- Confirm that the lot number of the control matches the lot number on the assay sheet.
- When using the MEK-9100 or MEK-9200 automated hematology analyzer, use the handy barcode reader provided with the hematology analyzer to read the 2D barcodes on the assay sheet of a Nihon Kohden genuine control. The assay values are entered into the hematology analyzer.
- Assay values of WBC 5 part differentials are for the case of optical adjustment with MEK-CAL HEMATOLOGY CALIBRATOR.

Measurement of a control is recommended for statistical quality control of the measurement accuracy of the analyzer. Also, it is recommended to use the average and upper and lower limits set by each laboratory to measure the control.

Frequency of Performing Quality Control

Determine according to the procedures established by each laboratory.

Control Used


MEK-5DL, MEK-5DN, MEK-5DH
HEMATOLOGY CONTROL for WBC 5 part differential

Measurement Principle

Refer to the analyzer operator's manual.

Procedure

1. Take the control out of the pharmaceutical refrigerator in which it is stored.
2. Allow the control to warm to room temperature (15 - 30°C, 59 - 86°F) for 15 minutes before mixing.

NOTE: Do not mix the control before it warms to room temperature. Doing so may cause the contents of the control to ooze out the top of the cap.

3. Mix the control.


NOTE

- Do not use a mechanical mixer.
- Do not mix the control roughly. If it bubbles, hemolysis is caused.

- 1) Hold the tube containing the control between the palms of the hands and roll it back and forth. Occasionally invert the tube so that the control is mixed thoroughly.

Roll the tube back and forth

Invert the tube

- 2) Continue to mix in the manner described in step 1) until the red cells are completely suspended. Tubes stored for a long time may require extra mixing.

4. Gently invert the tube until there is no clot of blood on the bottom before measurement.

NOTE: If mixing is incomplete, measurement results may be inaccurate.

5. Measure the control as instructed in the analyzer operator's manual.

NOTE: If the cap was opened for measurement, use a lint-free paper towel to wipe off any residual liquid on the cap or rim of the tube after measurement.

6. Return the control to the pharmaceutical refrigerator.

NOTE: Do not leave the control out at room temperature for an extended period of time after measurement. Doing so may affect the results of subsequent measurements.

Technical Information

Metrological Traceability of Values Assigned to the Control

For each parameter, the internationally recognized reference method is:

- 1) WBC and RBC: ICSH 1988: The assignment of values to fresh blood used for calibrating automated blood cell counters. Clin Lab Hematol.1988; 10:203-212
- 2) HGB: CLSI H15-A3
- 3) HCT: CLSI H7-A3
- 4) PLT: International Council for Standardization in Hematology Expert Panel on Cytometry and International Society of Laboratory Hematology Task Force on Platelet Counting. Platelet counting by RBC/platelet ratio method. A reference method. Am Journal of Clinical Pathology 115:460-464 2001

Assay Values

- Refer to the assay values for the specific instrument model on the provided assay sheet.
- Assay values are determined on well-maintained, properly calibrated instruments using the conformed reagents. Reagent differences, maintenance, operation technique, and calibration may contribute to inter-laboratory variation.
- Assigned values are presented as a Mean and Range. The Mean is derived from replicate testing on instruments operated and maintained according to the manufacturer's instructions. The Range is an estimate of variation between laboratories and also takes into account inherent imprecision of the method and expected biological variability of the control material.
- Assay values on a new lot of control should be confirmed before the new lot is put into routine use. Test the new lot when the instrument in good working order and quality control results on the old lot are acceptable. The laboratory's recovered mean should be within the assay range.
- For greater control sensitivity each laboratory should establish its own mean and acceptable range and periodically reevaluate the mean. The laboratory range may include values outside of the assay range. The user may establish assay values not listed on the assay sheet, if the control is suitable for the method.

NOTE: RDW-CV is displayed as RDW on the MEK-7222 and MEK-8222 automatic hematology analyzer.

Limitations

The performance of the control is assured only if it is properly stored and used as described in this operator’s manual. Incomplete mixing of the tube prior to use invalidates both the sample withdrawn and any remaining material in the tube.

Composition

The control contains human erythrocytes, simulated leukocytes and mammalian platelets suspended in a plasma-like fluid with preservatives.

Sterilization Method

The control is not intended to be sterilized or kept in a sterile environment.

Mathematical Approach upon which the Calculation of the Analytical Result is Made

The following parameters are calculated as follows ¹:

NE% = NE / WBC

LY% = LY / WBC

MO% = MO / WBC

EO% = EO / WBC

BA% = BA / WBC

IG% = IG / WBC

MCV = HCT / RBC

MCH = HGB / RBC

MCHC = HGB / HCT

MPV = PCT / PLT ²

¹ Without unit correction

² PCT: Platelet Crit

NOTE: RBC indices are calculated by rounding off to the whole number or the nearest tenth.

Environmental Conditions

Storage and Transport Environment

Temperature: 2 to 8°C (36 to 46°F)

Usage Environment

Temperature: 15 to 30°C (59 to 86°F)

Expiration Date

The expiration date is shown on the tube label and package.

Shelf Life Date After Opening

14 days (within the expiration date)

NOTE: Record the opening date at first use.

Package and Catalog Number

Model	Qty	Catalog Number
MEK-5DN	MEK-5DN (3.0 to 3.4 mL) × 2	5DN02
	MEK-5DN (3.0 to 3.4 mL) × 6	5DN06
MEK-5DL, MEK-5DH	MEK-5DL (3.0 to 3.4 mL) × 6, MEK-5DH (3.0 to 3.4 mL) × 6	5DLH6
MEK-5DL, MEK-5DN, MEK-5DH	MEK-5DL (3.0 to 3.4 mL) × 2, MEK-5DN (3.0 to 3.4 mL) × 2, MEK-5DH (3.0 to 3.4 mL) × 2	5DM02
	MEK-5DL (3.0 to 3.4 mL) × 4, MEK-5DN (3.0 to 3.4 mL) × 4, MEK-5DH (3.0 to 3.4 mL) × 4	5DM04

Disposal

⚠ WARNING

Dispose of the control according to your local laws and your facility's guidelines (including incineration, melt treatment, sterilization, disinfection and request for waste disposal) for disposing of infectious medical waste. Otherwise, it may affect the environment. If there is a possibility that the control may have been contaminated with infection, it may cause infection.

When disposing of the control, such as when the expiration date is past, follow the instructions on the SDS of the control.

Revision History

Edition	Date	Details	Code Number
1st Edition	06 Jan 2003	Initial issue	0604-018444
8th Edition	19 Apr 2024	Corrections	0614-908007A
9th Edition	28 Jan 2025	Corrections	0614-908007B

- NOTE
- The code number of this manual was changed from 0604-018444E to 0614-908007 when the manual was updated from 6th Edition to 7th Edition.
 - Changes made in the most recent edition are indicated by a bar in the left margin of each page.

Note for users in the territory of the EEA and Switzerland:
Any serious incident that has occurred in relation to the device must be reported to the European Representative designated by the manufacturer and the Competent Authority of the Member State of the EEA and Switzerland in which the user and/or patient is established.

Copyright Notice

The entire contents of this manual are copyrighted by Nihon Kohden. All rights are reserved.

 Manufacturer
NIHON KOHDEN CORPORATION
1-31-4 Nishiochiai, Shinjuku-ku,
Tokyo 161-8560, Japan
Phone +81 3-5996-8041
<https://www.nihonkohden.com/>

 European Representative
NIHON KOHDEN EUROPE GmbH
Raiffeisenstrasse 10, 61191 Rosbach, Germany
Phone +49 6003-827-0 Fax +49 6003-827-599
UK Responsible Person
NIHON KOHDEN UK LTD.
Unit 3, Heyworth Business Park,
Old Portsmouth Road, Peasmarsh,
Guildford, Surrey, GU3 1AF, UK
Phone +44 1483-333-1328

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Isporučeni materijali

- MEK-5DL, MEK-5DN, MEK-5DH
HEMATOLOGY CONTROL for WBC 5 part differential
- Ispitni list ¹

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- 1) WBC a RBC: ICSH 1988: The assignment of values to fresh blood used for calibrating automated blood cell counters. Clin Lab Hematol.1988; 10:203-212
- 2) HGB: CLSI H15-A3
- 3) HCT: CLSI H7-A3
- 4) PLT: International Council for Standardization in Hematology Expert Panel on Cytometry and International Society of Laboratory Hematology Task Force on Platelet Counting. Platelet counting by RBC/platelet ratio method. A reference method. Am Journal of Clinical Pathology 115:460-464 2001

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- 1) WBC ja RBC: ICSH 1988: The assignment of values to fresh blood used for calibrating automated blood cell counters. Clin Lab Hematol.1988; 10:203-212
- 2) HGB: CLSI H15-A3
- 3) HCT: CLSI H7-A3
- 4) PLT: International Council for Standardization in Hematology Expert Panel on Cytometry and International Society of Laboratory Hematology Task Force on Platelet Counting. Platelet counting by RBC/platelet ratio method. A reference method. Am Journal of Clinical Pathology 115:460-464 2001

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- 1) WBC et RBC: ICSH 1988: The assignment of values to fresh blood used for calibrating automated blood cell counters. Clin Lab Hematol.1988; 10:203-212
- 2) HGB : CLSI H15-A3
- 3) HCT : CLSI H7-A3
- 4) PLT : International Council for Standardization in Hematology Expert Panel on Cytometry and International Society of Laboratory Hematology Task Force on Platelet Counting. Platelet counting by RBC/platelet ratio method. A reference method. Am Journal of Clinical Pathology 115:460-464 2001

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- 1) WBC und RBC: ICSH 1988: The assignment of values to fresh blood used for calibrating automated blood cell counters. Clin Lab Hematol. 1988; 10:203-212
- 2) HGB: CLSI H15-A3
- 3) HCT: CLSI H7-A3
- 4) PLT: International Council for Standardization in Hematology Expert Panel on Cytometry and International Society of Laboratory Hematology Task Force on Platelet Counting. Platelet counting by RBC/platelet ratio method. A reference method. Am Journal of Clinical Pathology 115:460-464 2001

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- 1) WBC και RBC: ICSH 1988: The assignment of values to fresh blood used for calibrating automated blood cell counters. Clin Lab Hematol.1988; 10:203-212
- 2) HGB: CLSI H15-A3
- 3) HCT: CLSI H7-A3
- 4) PLT: International Council for Standardization in Hematology Expert Panel on Cytometry and International Society of Laboratory Hematology Task Force on Platelet Counting. Platelet counting by RBC/platelet ratio method. A reference method. Am Journal of Clinical Pathology 115:460-464 2001

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- 1) WBC un RBC: ICSH 1988: The assignment of values to fresh blood used for calibrating automated blood cell counters. Clin Lab Hematol.1988; 10:203-212
- 2) HGB: CLSI H15-A3
- 3) HCT: CLSI H7-A3
- 4) PLT: International Council for Standardization in Hematology Expert Panel on Cytometry and International Society of Laboratory Hematology Task Force on Platelet Counting. Platelet counting by RBC/platelet ratio method. A reference method. Am Journal of Clinical Pathology 115:460-464 2001

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- 1) WBC și RBC: ICSH 1988: The assignment of values to fresh blood used for calibrating automated blood cell counters. Clin Lab Hematol.1988; 10:203-212
- 2) HGB: CLSI H15-A3
- 3) HCT: CLSI H7-A3
- 4) PLT: International Council for Standardization in Hematology Expert Panel on Cytometry and International Society of Laboratory Hematology Task Force on Platelet Counting. Platelet counting by RBC/platelet ratio method. A reference method. Am Journal of Clinical Pathology 115:460-464 2001

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- 1) WBC и RBC: ICSH 1988: The assignment of values to fresh blood used for calibrating automated blood cell counters. Clin Lab Hematol.1988; 10:203-212
- 2) HGB: CLSI H15-A3
- 3) HCT: CLSI H7-A3
- 4) PLT: International Council for Standardization in Hematology Expert Panel on Cytometry and International Society of Laboratory Hematology Task Force on Platelet Counting. Platelet counting by RBC/platelet ratio method. A reference method. Am Journal of Clinical Pathology 115:460-464 2001

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- 1) WBC in RBC: ICSH 1988: The assignment of values to fresh blood used for calibrating automated blood cell counters. Clin Lab Hematol.1988; 10:203-212
- 2) HGB: CLSI H15-A3
- 3) HCT: CLSI H7-A3
- 4) PLT: International Council for Standardization in Hematology Expert Panel on Cytometry and International Society of Laboratory Hematology Task Force on Platelet Counting. Platelet counting by RBC/platelet ratio method. A reference method. Am Journal of Clinical Pathology 115:460-464 2001

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- 1) WBC y RBC: ICSH 1988: The assignment of values to fresh blood used for calibrating automated blood cell counters. Clin Lab Hematol.1988; 10:203-212
- 2) HGB: CLSI H15-A3
- 3) HCT: CLSI H7-A3
- 4) PLT: International Council for Standardization in Hematology Expert Panel on Cytometry and International Society of Laboratory Hematology Task Force on Platelet Counting. Platelet counting by RBC/platelet ratio method. A reference method. Am Journal of Clinical Pathology 115:460-464 2001

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