



Paste for EEG

Z-101N

General

The Z-101N Paste for EEG is a conductive paste used with surface electrodes to reduce the contact impedance between electrodes and the patient's skin when performing EEG measurement.

Safety Information

⚠ CAUTION A caution alerts the user to possible injury or problems with the instrument associated with its use or misuse such as instrument malfunction, instrument failure, damage to the instrument, or damage to other property.

Pay attention to all safety information in the Operator's Manual or Installation Guide.

⚠ CAUTION

- Use the paste only for EEG measurement.
- If the paste contacts the eyes or mouth, wash thoroughly and immediately with water and see a physician.
- If rash, redness, itch, swelling or abnormality appears on skin, wash the paste off immediately and do not use it.
- Do not leave the paste near the patient or in reach of children.
- If the paste is swallowed, see a physician.

⚠ CAUTION

Dispose of Nihon Kohden products according to your local laws and your facility's guidelines for waste disposal. Otherwise, it may affect the environment. If there is a possibility that the product may have been contaminated with infection, dispose of it as medical waste according to your local laws and your facility's guidelines for medical waste. Otherwise, it may cause infection.

This Safety and Performance Information is an extract from the general and safety information sections of the most recent edition of Operator's Manual or Installation Guide. Therefore, the contents of your Operator's Manual or Installation Guide may differ from those of this Safety and Performance Information. For detailed operating procedures, follow the instructions of your Operator's Manual or Installation Guide.



Manufacturer

NIHON KOHDEN CORPORATION

1-31-4 Nishi-chiai, Shinjuku-ku,

Tokyo 161-8560, Japan

Phone +81 3-5996-8041

<https://www.nihonkohden.com/>



European Representative

NIHON KOHDEN EUROPE GmbH

Raiffeisenstrasse 10, 61191 Rosbach, Germany

Phone +49 6003-827-0 Fax +49 6003-827-599

1st Edition: 12 Mar 2019

2nd Edition: 13 Sep 2024

