Celltac chemi Rinse solution RN-410V

English

General

Intended Purpose

For in vitro diagnostic use only.

The Celltac chemi Rinse solution is used for in vitro diagnosis as a rinse solution for filling the internal flow path of the Nihon Kohden Clinical Chemistry Analyzer, maintaining the pressure of the syringe piston, and rinsing blood and reagents. This device is intended for laboratory professional use only.

Read the analyzer operator's manual together with this manual before and during use.

- NOTE**
- Use the Celltac chemi Rinse solution with the specified analyzers only.
 - Read the SDS (Safety Data Sheet) carefully before use. The SDS is available from your Nihon Kohden representative.

Devices Intended for Use in Combination with
CHM-4100 clinical chemistry analyzer

Materials Provided and Materials Required

Materials Provided

RN-410V Celltac chemi Rinse solution

Materials Required (Not Provided with the Celltac chemi Rinse solution)




None

Intended Users

For laboratory professional use only, in laboratories with suitable equipment for clinical chemistry testing. Qualified personnel, e.g. laboratory technicians trained in clinical chemistry analysis techniques, will be able to use according to the Operator's Manual.

Symbols

The following symbols are used with the Celltac chemi Rinse solution. The descriptions of each symbol are given in the table below.

Symbol	Description	Symbol	Description
	Use by		Operator's manual; operating instructions
	Lot number		In vitro diagnostic medical device
	Catalogue number		Exclamation mark ¹
	Fragile		Unique device identifier
	Keep away from sunlight		Model number
	Temperature limits		Manufacturer
	This way up		Authorized representative
	Keep away from rain		The CE mark is a protected conformity mark of the European Union.
	Stacking limit by number ("n" is the limiting number)		Importer
	Caution		

¹ The label is affixed in accordance with GHS requirements. For details, refer to "Hazards Identification" (p. 2).

Safety Information

	CAUTION	A caution alerts the user to possible injury or problems with the instrument associated with its use or misuse such as instrument malfunction, instrument failure, damage to the instrument, or damage to other property.
---	----------------	---

Pay attention to all safety information in this operator's manual.

	CAUTION
<ul style="list-style-type: none"> • Do not swallow the Celltac chemi Rinse solution. If swallowed, rinse the mouth immediately. Do not force vomiting. See a physician. • If the Celltac chemi Rinse solution contacts the eyes or mouth, wash thoroughly and immediately with plenty of water and see a physician. • If the Celltac chemi Rinse solution contacts the skin, wash with plenty of water. 	

NOTE: Do not use the Celltac chemi Rinse solution if it is past the expiration date on the package or if the Celltac chemi Rinse solution is stored under unspecified conditions.

Hazards Identification

Signal Word

Warning

Hazard Statement

- H317 May cause an allergic skin reaction.
- H412 Harmful to aquatic life with long lasting effects

Precautionary Statement Prevention

- P261 Avoid breathing dust/fume/gas/mist/vapours/spray.
- P272 Contaminated work clothing should not be allowed out of the workplace.
- P273 Avoid release to the environment.
- P280 Wear protective gloves/protective clothing/eye protection/face protection.

Precautionary Statement Response

- P302+P352
IF ON SKIN: Wash with plenty of water.
- P333+P313
If skin irritation or rash occurs: Get medical advice/attention.
- P362+P364
Take off contaminated clothing and wash it before reuse.

Precautionary Statement Disposal

- P501 Dispose of contents/container in accordance with local and national regulations.

Reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1) $\geq 0.0015\%$

Using the Celltac chemi Rinse solution

Refer to the analyzer operator's manual for details.

- NOTE
- Use the Celltac chemi Rinse solution without dilution or mixing.
 - Use the Celltac chemi Rinse solution when its temperature is between 15 and 30°C (59 and 86°F).
 - Do not add the Celltac chemi Rinse solution to another rinse solution bottle.

Technical Information

Composition

0.003% or less of Reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1) is contained as preservative.

Sterilization Method

The Celltac chemi Rinse solution is not intended to be sterilized or kept in a sterile environment.

Environmental Conditions

Storage and Transport Environment

Temperature: 1 to 30°C (34 to 86°F)

Usage Environment

Temperature: 15 to 30°C (59 to 86°F)

Expiration Date

The expiration date is shown on the container label and package.

Shelf Life Date After Opening

Use the Celltac chemi Rinse solution within 90 days of opening.

Package and Catalog Number

Model	Qty	Catalog Number
RN-410V	250 mL × 1	T561
	250 mL × 3	T561A

Disposal

⚠ CAUTION

Dispose of the Celltac chemi Rinse solution according to your local laws and your facility's guidelines (including incineration, melt treatment, sterilization, disinfection and request for waste disposal) for disposing of infectious medical waste. Otherwise, it may affect the environment. If there is a possibility that the Celltac chemi Rinse solution may have been contaminated with infection, it may cause infection.

When disposing of the Celltac chemi Rinse solution, such as when the expiration date is past, follow the instructions on the SDS of the Celltac chemi Rinse solution.

Revision History

Edition	Date	Details	Code Number
1st Edition	05 Aug 2016	Initial issue	0614-907557
4th Edition	27 Apr 2023	Symbols added	0614-907557D
5th Edition	06 Mar 2026	Symbols and notes are added	0614-907557E

NOTE: Changes made in the most recent edition are indicated by a bar in the left margin of each page.

Note for users in the territory of the EEA and Switzerland:
Any serious incident that has occurred in relation to the device must be reported to the European Representative designated by the manufacturer and the Competent Authority of the Member State of the EEA and Switzerland in which the user and/or patient is established.

Copyright Notice

The entire contents of this manual are copyrighted by Nihon Kohden. All rights are reserved.

 Manufacturer
NIHON KOHDEN CORPORATION
1-31-4 Nishiochiai, Shinjuku-ku,
Tokyo 161-8560, Japan
Phone +81 3-5996-8041
<https://www.nihonkohden.com/>

 European Representative
NIHON KOHDEN EUROPE GmbH
Raiffeisenstrasse 10, 61191 Rosbach, Germany
Phone +49 6003-827-0 Fax +49 6003-827-599
UK Responsible Person
NIHON KOHDEN UK LTD.
Unit 3, Heyworth Business Park,
Old Portsmouth Road, Peasmarsh,
Guildford, Surrey, GU3 1AF, UK
Phone +44 14-8333-1328

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