# Celltac chemi HbA1c N

## HA-420W, HA-421W

English

## General

### Intended Purpose

*For in vitro diagnostic use only.*

Celltac chemi HbA1c N is for in vitro diagnostic use in the automated quantitative measurement of glycosylated hemoglobin A1c (HbA1c, Hemoglobin A1c) in EDTA or sodium fluoride treated whole blood, as an aid to diagnosis in the detection and evaluation of diabetes mellitus. For use with CHM-4100 clinical chemistry analyzer (only when software Ver.03-01 or later is installed) and MEK-1303 automated hematology and clinical chemistry analyzer. The measurement of HbA1c is recommended as a marker of long-term metabolic control in persons with Diabetes Mellitus. The Celltac chemi HbA1c N is intended for laboratory professional use only. Testing population are general population patients with suspected hyperglycemia or monitored for diabetes.

Read the analyzer operator's manual together with this manual before and during use.

- NOTE**
- Nihon Kohden representative will replace a product which proves to be defective until the expiration date shown on the package, provided that the product is used as prescribed by the operating instructions in the operator's manuals.
  - Overall judgement must be performed by the physician, referring to the analysis result, clinical findings and other examination results.
  - Use a test cartridge for a single test. Do not reuse it.
  - Use the product with the specified analyzers only.
  - Read the SDS (Safety Data Sheet) carefully before use. The SDS is available from your Nihon Kohden representative.
  - The SSP (Summary of Safety and Performance) for the Celltac chemi HbA1c N can be obtained from the following URL:  
<https://ec.europa.eu/tools/eudamed>

### Devices Intended for Use in Combination with

- CHM-4100 clinical chemistry analyzer (only when software Ver. 03-01 or later is installed)

- MEK-1303 automated hematology and clinical chemistry analyzer

## Materials Provided and Materials Required

### Materials Provided

Celltac chemi HbA1c N contains the following reagents used for HbA1c measurement in a single cartridge.

- A1c diluent
- A1c latex reagent
- A1c antibody reagent
- A1c rinse solution

### Materials Required (Not Provided with the Reagent)

- YZ-005B1 HA-410V CAL
- HbA1c control (commercially available control materials)

## Summary and Explanation

Measurement of hemoglobin A1c (HbA1c) is accepted as a method to measure long-term glucose control in patients with diabetes mellitus. Measurement of HbA1c is an important tool for monitoring the efficacy of dietary control and therapy during treatment of diabetes mellitus. Long term treatment of the disease emphasizes control of blood glucose levels in preventing the acute complications of ketosis and hyperglycemia. In addition, long term complications such as retinopathy, neuropathy, and cardiovascular disease can be minimized if blood glucose levels are effectively controlled. The process of conversion from hemoglobin A to hemoglobin A1c depends on the blood glucose concentration. Since the average life of a red blood cell is 120 days, measurement of haemoglobin A1c can reflect the mean daily blood glucose concentration over the preceding two to three months and provide a much better indication of glycemic control than blood or urinary glucose measurement.

## Analyte or Marker

Celltac chemi HbA1c N is intended to measure the ratio (%) of glycohemoglobin A1c (HbA1c) in total hemoglobin (Hb).

## Target Treated Population

Testing population are patients with suspected hyperglycemia or monitored for diabetes. General population (adult and pediatric).

## Specimen Collection and Preparation


- Use a sample of whole blood.
- When collecting whole blood specimen, use the K<sub>2</sub>-EDTA, K<sub>3</sub>-EDTA, Na<sub>2</sub>-EDTA or sodium fluoride blood collection tube.
- Use the sample within the following period after drawing blood.
  - When storing the sample at 4°C (39°F) or less: 5 days
  - When storing the sample at room temperature (around 25°C [77°F]): 2 days)
- For specimen collection and transportation, follow the instructions for use of the EDTA-blood collection tube or sodium fluoride blood collection tube and general recommendations for blood sampling.



## Intended Users

This device is intended to be for professional use only. Qualified personnel, e.g. laboratory technicians trained in hematology analysis techniques will be able to use according to the operator's manual.

## Symbols

The following symbols are used with the measurement kit. The descriptions of each symbol are given in the table below.

Symbol	Description
	Use by
	Lot number
	Catalogue number
	Do not reuse
	Fragile
	Keep away from sunlight
	Temperature limits
	Contains sufficient for <n> tests
	Caution
	Operator's manual; operating instructions
	Model number

Symbol	Description
	Unique device identifier
	In vitro diagnostic medical device
	Exclamation mark <sup>1</sup>
	Manufacturer
	Authorized representative
	The CE mark is a protected conformity mark of the European Union. The four digits after the CE mark indicate the identification number of the Notified Body involved in assessing the product's conformity as a medical device.
	Importer

<sup>1</sup> The label is affixed in accordance with GHS requirements. For details, refer to "Hazards Identification" (p. 3).

## Safety Information

- ⚠ WARNING** A warning alerts the user to the possible injury or death associated with the use or misuse of the instrument.
- ⚠ CAUTION** A caution alerts the user to possible injury or problems with the instrument associated with its use or misuse such as instrument malfunction, instrument failure, damage to the instrument, or damage to other property.

Pay attention to all safety information in this operator's manual.

- ⚠ WARNING**
- Always wear rubber gloves to protect yourself from infection when handling and measuring blood samples and test cartridge. The waste fluid may leak from the test cartridge after use.

- ⚠ CAUTION**
- Do not swallow the reagent. If swallowed, rinse the mouth immediately. Do not force vomiting. See a physician.
  - If the reagent contacts the eyes, mouth or skin, wash immediately with plenty of water and see a physician.

- NOTE**
- Do not use the product if it is past the expiration date on the package or if the product is stored under unspecified conditions.
  - Store the product inside a polystyrene container and put it in a pharmaceutical refrigerator (2 to 8°C, 36 to 46°F) which allows temperature adjustment. Do not store the product in a domestic refrigerator. This may freeze the reagent.
  - Do not freeze the product. Do not use it if it is frozen.
  - The usage environment of the product is between 15 and 30°C (59 and 86°F).
  - Be careful when handling the reagent because it contains a small amount of toxic sodium azide as a preservative.
  - The outer film of this product is a gas barrier film. After removing the outer film, do not store the product in an environment where substances that release carbon dioxide gas (e.g. dry ice) are present. The carbon dioxide gas may affect the product, resulting in inaccurate measurements and improper diagnoses.
  - Before removing the outer film of the product, check that the film is not damaged. If the film is damaged, do not use the product.

## Hazards Identification

### Signal Word

Warning

### Hazard Statement

- H317 May cause an allergic skin reaction.
- EUH032 Contact with acids liberates very toxic gas.
- EUH071 Corrosive to the respiratory tract.

### Precautionary Statement Prevention

- P261 Avoid breathing dust/fume/gas/mist/vapours/spray.
- P272 Contaminated work clothing should not be allowed out of the workplace.
- P280 Wear protective gloves/protective clothing/eye protection/face protection.

### Precautionary Statement Response

- P302 + P352 IF ON SKIN: Wash with plenty of water.
- P333 + P313 If skin irritation or rash occurs: Get medical advice/attention.
- P362 + P364 Take off contaminated clothing and wash it before reuse.

### Precautionary Statement Disposal

- P501 Dispose of contents/container in accordance with local and national regulations.

Reaction mass of 5-Chloro-2-methyl-4-isothiazolin-3-one and 2-Methyl-2H -isothiazol-3-one (3:1)  $\geq 0.0015\%$

2-methylisothiazol-3(2H)-one  $\geq 0.0015\%$

Sodium azide  $< 0.1\%$

## Measurement

NOTE: Refer to the analyzer operator's manual for details on measurement.


### Measurement Principle

This product uses the latex agglutination immunoassay method to measure HbA1c concentration as a percentage of total hemoglobin (Hb and HbA1c) in whole blood. Whole blood samples are hemolyzed by A1c diluent. Total hemoglobin in the hemolyzed blood binds with the latex particles. The amount of binding is proportional to the relative concentration of both Hb and HbA1c in the blood. The A1c antibody reagent binds to the particle bound HbA1c and agglutination takes place. The measured absorbance is proportional to the HbA1c which is bound to particles, which in turn is proportional to the percentage of HbA1c in the sample.

### Procedure

1. Prepare the analyzer. Refer to the analyzer operator's manual.
2. Take a test cartridge out of the refrigerator. Wipe off condensation with a dry cloth.
3. Set the test cartridge into the analyzer. Gently invert the test several times. If the agitation is not enough, the result may be inaccurate.

Invert the cartridge more than 3 times.

The analyzer starts reading the QR code on the cartridge, which includes the calibration curve information of the reagent. If an error message appears, refer to the analyzer operator's manual.

4. Measure a sample of whole blood.

Use a sample with a hemoglobin concentration of 7 to 20 g/dL. Otherwise, the test result may be inaccurate.

When using the whole blood sample, gently invert and mix the sample before use because the red cell components settle easily and become unhomogenous.

The measurement results are reported in %NGSP. This unit corresponds to HbA1c concentration as a percentage of total hemoglobin. The following equation shows the relationship between NGSP and IFCC.


NGSP value (%) =  $0.0915 \times \text{IFCC value (mmol/mol)} + 2.15$

- NOTE
- No dirt or damage on the test cartridge
  - No leakage from the test cartridge

- Before setting the test cartridge, check that the reagent is at the bottom of the wells and that there are no bubbles at the bottom of the wells. If the reagent is at the top of the wells or there are bubbles at the bottom of the wells, gently shake the cartridge to make the reagent go down. Otherwise the test result may be inaccurate.

- Do not use a test cartridge which fell on the floor because the light path surface may be dirty.
- Do not scratch the test cartridge.
- Do not touch the QR code and optical path side of the cartridge.

- The correlation between measurements obtained using this product and measurements obtained using existing products has been statistically verified. However variations may be observed in practice as a result of differences in measurement methods, facilities and equipment.
- False reactions have been observed in cases of patients with hemoglobinopathy or high gamma globulin levels. If the measurement value is too high or too low for the patient condition or the measurement result is suspicious, remeasure the sample or measure the sample by another measurement method.

## Calibration

NOTE: Contact your Niohn Kohden representative when HbA1c calibration is required.

Perform calibration in the following cases.

- Calibration is required as a result of regular quality control.
- Calibration is required as a result of quality control performed after maintenance.
- Other cases where calibration is required.

## Reagents Used

Use an HA-410V CAL in combination with the HA-420W or HA-421W Celltac chemi HbA1c N.

## Frequency of Performing Calibration

Determine according to inspections performed by each laboratory.

## Procedure

Refer to the analyzer operator's manual or service manual.

## Quality Control

Measurement of a control is recommended for statistical quality control of the measurement accuracy of the analyzer. Also, it is recommended to use the average and upper and lower limits set by each laboratory to measure the control.

## Frequency of Performing Quality Control

Determine according to the procedures established by each laboratory.

## Control Used

Commercially available control materials.

## Procedure

Refer to the analyzer operator's manual.

## Variable Factors

Variable factors that should be considered when setting appropriate upper and lower limits for quality control include the following.

- Calibrator lots
- Control lots, etc.
- Reagent lots
- Analyzer used for the control measurements
- Measurement frequency

## Technical Information

### Metrological Traceability of Values Assigned to the Calibrator

Assigned values of HbA1c traceable to an internationally recognized reference material: JCCRM-411

The measurement results are reported in %NGSP. This unit corresponds to HbA1c concentration as a percentage of total hemoglobin. The following equation shows the relationship between NGSP and IFCC.

$$\text{NGSP value (\%)} = 0.0915 \times \text{IFCC value (mmol/ mol)} + 2.15$$

Results are traceable to HbA1c calibrator for NGSP certification (JCCAL)/Certified Reference Material for Measurement of HbA1c (JCCRM).

## Quality Control, Analytical Performance Characteristics

### Sensitivity Test

When the standard solution containing about 9.5% of HbA1c concentration is used as the sample solution, the absorbance change is within the range of 0.16 to 0.57.

### Accuracy Test

When the control sample of known concentration is measured, the measured value is within  $\pm 10\%$  of the known concentration.

### Repeatability Test

When the control samples of known concentration is measured, the CV of the measured values is 5% or less.

### Measurement Range (Measuring Interval)

The lower limit of the measuring interval is defined by the limit of quantitation (LoQ). The upper limit of the measuring interval was established through a linearity study.

HbA1c: 4 to 13 % (NGSP)<sup>1</sup>

<sup>1</sup> Confirmed on the MEK-1303 and CHM-4100

### Expected Values

- Normal range:  $< 5.8\%$  <sup>1, 2</sup>
- Prediabetes: 5.8 to 6.4% <sup>1, 2</sup>
- Diabetes:  $\geq 6.5\%$  <sup>1, 2</sup>
- Cut-off Value: 6.5% <sup>3</sup>

<sup>1</sup> HbA1c: A Review of Analytical and Clinical Aspects. Ann lab Med 2013; 33: 393-400

<sup>2</sup> American Diabetes Association: Standards of Medical Care in Diabetes—2010. Diabetes Care 2010; 33 (Suppl. 1): S11-S61

<sup>3</sup> Journal of Diabetes Investigation Volume 3 Issue 1 February 2012, WHO Report WHO/NMH/CHP/CPM/11.1: Use of Glycated Hemoglobin (HbA1c) in the Diagnosis of Diabetes Mellitus & Diabetes Care Jan 2018, Volume 41, Supplement 1.

**NOTE:** Do not make decisions about drug selection or application of hypoglycemic drugs based solely on results obtained using this product. Instead, make an overall decision based on the medical guidelines of your facility considering a range of factors such as blood sugar levels and other test results.

### Batch Variation

Three lot of the Celltac Chemi HbA1c N (HA-420W/421W) on the MEK-1303 and CHM-4100 was used to evaluate between lot. Two levels of Liquichek Diabetes Control was 5 repeated for 3 runs. The %CVs of the measurements were calculated.

Model	Control Level 1	Control Level 2
MEK-1303	1.5%	1.9%
CHM-4100	0.0%	0.4%

## Composition

Test cartridge composition:

- A1c diluent: Colorless transparent liquid (preservative agent: Proclin300, sodium azide)
- A1c latex reagent: Milky white liquid (buffering agent, polystyrene latex, preservative agent: Proclin950)
- A1c antibody reagent: Colorless transparent liquid (buffering agent, anti-human HbA1c mouse monoclonal antibody: 1.6 $\pm$ 0.5  $\mu$ g per cartridge, preservative agent: Proclin950)
- A1c rinse solution: Colorless transparent liquid (preservative agent: Proclin300, sodium azide)

## Sterilization Method

The HA-420W or HA-421W Celltac chemi HbA1c N is not intended to be sterilized or kept in a sterile environment.

## Interfering Substances or Limitations

### Interfering Substances

The following interfering substances have been confirmed to have no effect on test results below the listed concentrations.

Substance	Max. Concentration
Rheumatoid factor	1,062 IU/mL
Bilirubin C (conjugated type)	34.9 mg/dL
Bilirubin F (free type)	31.8 mg/dL
Chyle	2,769 FTU

### Cross-reactivity

The following hemoglobin derivatives resulting from the following concentrations of the following substances will not affect the test result.

- Unstable hemoglobin from up to 2,000 mg/dL glucose
- Carbamylated hemoglobin from up to 80 mg/dL sodium cyanate
- Acetylated hemoglobin from up to 150 mg/dL acetylsalicylic acid
- Acetaldehyde-hemoglobin from up to 100 mg/dL acetaldehyde
- Sodium nitrite induced methemoglobin at a concentration of 60% or less

## Analytical Performance Characteristics

Values shown in this section are representative performance data of the Celltac chemi HbA1c N and may differ from values acquired in individual facilities.

### Analytical Sensitivity

To determine Limit of Detection (LoD), Limit of Blank (LoB), and Limit of Quantification (LoQ) of the Celltac chemi HbA1c N on the MEK-1303 and CHM-4100.

Reference guidance: CLSI EP17-A2

Model	LoB	LoD	LoQ
MEK-1303	1.5%	1.6%	1.9%
CHM-4100	1.4%	1.5%	2.0%

### Precision

One lot of Celltac Chemi HbA1c N (HA-420W/421W) was used to evaluate precision. Two levels of Liquichek Diabetes Control (Bio-Rad Laboratories, Inc.) was repeated twice in the morning and evening for 20 days.

Reference guidance: CLSI EP05-A3

Model	Sample	Mean	n	Repeatability		Within-laboratory precision	
				SD	%CV	SD	%CV
MEK-1303	Control Level 1	5.34%	80	0.052	1.0%	0.121	2.3%
	Control Level 2	9.30%	80	0.054	0.6%	0.109	1.2%
CHM-4100	Control Level 1	5.19%	80	0.045	0.9%	0.089	1.7%
	Control Level 2	9.27%	80	0.064	0.7%	0.084	0.9%

### Accuracy, Trueness

Primary Practical Sample Certified Reference Material for HbA1c Measurement JCCRM411-3 levels 1 to 5 were measured in triplicate across 3 different reagent lots on the MEK-1303 and CHM-4100 analyzers to confirm that each reagent batch was measuring HbA1c accurately.

Model	Reference Material	Certified value(%)	Reagent Lot		
			1	2	3
MEK-1303	JCCRM Level 1	5.10	98.7%	102.0%	102.0%
	JCCRM Level 2	5.77	100.5%	102.3%	102.3%
	JCCRM Level 3	7.39	100.1%	102.4%	101.5%
	JCCRM Level 4	9.60	99.3%	102.8%	102.1%
	JCCRM Level 5	11.98	99.1%	101.8%	100.7%
CHM-4100	JCCRM Level 1	5.10	100.0%	102.0%	101.3%
	JCCRM Level 2	5.77	100.5%	101.1%	100.5%
	JCCRM Level 3	7.39	101.5%	101.5%	101.0%
	JCCRM Level 4	9.60	100.7%	101.7%	102.4%
	JCCRM Level 5	11.98	105.5%	104.1%	104.9%

### Linearity

Samples were measured in quadruplicate using the MEK-1303 and CHM-4100 analysers to confirm linearity across the measuring range.

Reference guidance: CLSI EP06-A

This kit measuring interval obtained was below:

Model	Linearity Measuring Interval
MEK-1303	3.50 to 13.25%
CHM-4100	3.20 to 13.50%

### Correlation

To determine the correlation between the Celltac chemi HbA1c N on the MEK-1303 and Celltac chemi HbA1c N on the CHM-4100 and HPLC

Correlation coefficient and slope are calculated by least squares method

Reference guidance: CLSI EP09-A3

Model	n	Correlation Coefficient	Slope
MEK-1303 v.s. HPLC(Other company)	116	0.995	0.9581
CHM-4100 v.s. HPLC(Other company)	116	0.987	0.9448

- CLSI. *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. CLSI document EP17-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2012.
- CLSI. *Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline—Third Edition*. CLSI document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
- CLSI. *Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline*. CLSI document EP06-A. Wayne, PA: Clinical and Laboratory Standards Institute; 2003.
- CLSI. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. CLSI document EP09-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.

## Mathematical Approach upon which the Calculation of the Analytical Result is Made

Calibration information is determined at the factory for each reagent lot. The calibration information is recorded in the QR code, and the analyzer reads the information before measurement and calculates the HbA1c value.

## Environmental Conditions

### Storage and Transport Environment

Temperature: 2 to 8°C (36 to 46°F)

## Usage Environment


Temperature: 15 to 30°C (59 to 86°F)

## Expiration Date


The expiration date is shown on the package.

Lifetime: 18 months


## Package and Catalog Number

HA-420W

HA-421W

Test cartridge

Model	Qty	Catalog Number
HA-420W	50 (25 test cartridges × two boxes)	HA-420W
HA-421W	25	HA-421W

## Disposal

### WARNING

Dispose of the test cartridge according to your local laws and your facility's guidelines (including incineration, melt treatment, sterilization, disinfection and request for waste disposal) for **disposing of infectious medical waste**. Otherwise, it may affect the environment. If there is a possibility that the product may have been contaminated with infection, it may cause infection.

When disposing of the test cartridges, such as when the expiration date is past, follow the instructions on the SDS of the Celltac chemi HbA1c N.

## Revision History

Edition	Date	Details	Code Number
1st Edition	31 Oct 2018	Initial issue	0654-905277
6th Edition	08 Jan 2025	Corrections	0614-907957C
7th Edition	06 Mar 2026	Symbols and notes are added	0614-907957D

- NOTE
- The code number of this manual was changed from 0654-905277B to 0614-907957A when the manual was updated from 3rd Edition to 4th Edition.
  - Changes made in the most recent edition are indicated by a bar in the left margin of each page.

### Trademark

“QR Code” is registered trademark of DENSO WAVE INCORPORATED.

Note for users in the territory of the EEA and Switzerland:  
Any serious incident that has occurred in relation to the device must be reported to the European Representative designated by the manufacturer and the Competent Authority of the Member State of the EEA and Switzerland in which the user and/or patient is established.

### Copyright Notice

The entire contents of this manual are copyrighted by Nihon Kohden. All rights are reserved.

 Manufacturer  
**NIHON KOHDEN CORPORATION**  
1-31-4 Nishiochiai, Shinjuku-ku,  
Tokyo 161-8560, Japan  
Phone +81 3-5906-8041  
<https://www.nihonkohden.com/>

**EU REP** European Representative  
**NIHON KOHDEN EUROPE GmbH**  
Raiffeisenstrasse 10, 61191 Rosbach, Germany  
Phone +49 6003-827-0 Fax +49 6003-827-599  
UK Responsible Person  
**NIHON KOHDEN UK LTD.**  
Unit 3, Heyworth Business Park,  
Old Portsmouth Road, Peasmarsh,  
Guildford, Surrey, GU3 1AF, UK  
Phone +44 14-8333-1328

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- CLSI. *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. CLSI-Dokument EP17-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2012.
- CLSI. *Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline—Third Edition*. CLSI-Dokument EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
- CLSI. *Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline*. CLSI document EP06-A. Wayne, PA: Clinical and Laboratory Standards Institute; 2003.
- CLSI. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. CLSI-Dokument EP09-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.

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- CLSI. *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. CLSI document EP17-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2012.

- CLSI. *Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline—Third Edition*. CLSI document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
- CLSI. *Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline*. CLSI document EP06-A. Wayne, PA: Clinical and Laboratory Standards Institute; 2003.
- CLSI. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. CLSI document EP09-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.

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- CLSI. *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. CLSI document EP17-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2012.
- CLSI. *Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline—Third Edition*. CLSI document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.

- CLSI. *Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline*. CLSI document EP06-A. Wayne, PA: Clinical and Laboratory Standards Institute; 2003.
- CLSI. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. CLSI document EP09-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.

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- CLSI. *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. CLSI document EP17-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2012.
- CLSI. *Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline—Third Edition*. CLSI document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
- CLSI. *Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline*. CLSI document EP06-A. Wayne, PA: Clinical and Laboratory Standards Institute; 2003.

- CLSI. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. CLSI document EP09-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.

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- CLSI. *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. Documento EP17-A2 del CLSI. Wayne, PA: Clinical and Laboratory Standards Institute; 2012.
- CLSI. *Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline—Third Edition*. Documento EP05-A3 del CLSI. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.

- CLSI. *Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline*. Documento EP06-A del CLSI. Wayne, PA: Clinical and Laboratory Standards Institute; 2003.
- CLSI. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. Documento EP09-A3 del CLSI. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.

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