HEMATOLOGY CALIBRATOR MEK-CAL

English

General

Intended Purpose

For in vitro diagnostic use only.

HEMATOLOGY CALIBRATOR is for in vitro diagnostic use as a value assigned whole blood calibrator designed to calibrate the values of Nihon Kohden hematology analyzers. HEMATOLOGY CALIBRATOR is not used for diagnosis or assistance in diagnosis, prognosis, prediction of disease or patient's physiological condition. Measuring this HEMATOLOGY CALIBRATOR ensures that the device is functioning according to its intended use.

Calibrate parameters are below:

- 1) WBC: White Blood Cell Count
- 2) RBC: Red Blood Cell Count
- 3) HGB: Hemoglobin Concentration
- 4) HCT: Hematocrit Percent
- 5) MCV: Mean Corpuscular Volume
- 6) MCH: Mean Corpuscular Hemoglobin
- 7) MCHC: Mean Corpuscular Hemoglobin Concentration
- 8) PLT: Platelet Count
- 9) MPV: Mean Platelet Volume
- 10) RDW: Red Blood Cell Distribution Width

Read the analyzer operator's manual together with this manual before and during use. Also refer to the assay sheet available on the Nihon Kohden website (See the "Materials Provided" section).

- NOTE**
- Use the calibrator with the specified analyzers only.
 - Read the SDS (Safety Data Sheet) carefully before use. The SDS is available from your Nihon Kohden representative and Nihon Kohden Website (<https://www.nihonkohden.com>).

Devices Intended for Use in Combination with

MEK-1301, MEK-1302, MEK-1303, MEK-1305, MEK-6400, MEK-6410, MEK-6420, MEK-6450, MEK-6500, MEK-6510, MEK-6550, MEK-7222, MEK-7300, MEK-8222, MEK-9100, MEK-9200

0614-908008B

Materials Provided and Materials Required

Materials Provided

- MEK-CAL
HEMATOLOGY CALIBRATOR
- Assay sheet ¹

¹ The assay sheet can be found on the Nihon Kohden website (<https://www.nihonkohden.com/technology/assaysheet/assaysheet.html>). Print it on a laser printer with print quality set to 1200×1200 dpi or higher.

Materials Required (Not Provided with the Calibrator)



- Protective gloves
- Protective glasses
- Waste container for biohazardous materials

Intended Users

For laboratory professional use only, in laboratories with suitable equipment for hematological testing. Qualified personnel, e.g. laboratory technicians trained in hematology analysis techniques, will be able to use according to the Operator's Manual.

Symbols

The following symbols are used with the calibrator. The descriptions of each symbol are given in the table below.

Symbol	Description	Symbol	Description
	Caution		Biological risks
	Operator's manual; operating instructions		This way up
	Catalogue number		Calibrator
	Lot number		Open tube stability
	Use by		Unique Device Identifier
	In vitro diagnostic medical device		The CE mark is a protected conformity mark of the European Union. The four digits after the CE mark indicate the identification number of the Notified Body involved in assessing the product's conformity as a medical device.
	Manufacturer		
	Authorized representative in the European Community/ European Union		Model number

Safety Information

⚠ WARNING A warning alerts the user to the possible injury or death associated with the use or misuse of the instrument.

Pay attention to all safety information in this operator's manual.

⚠ WARNING

- **POTENTIALLY BIOHAZARDOUS MATERIAL.** For in vitro diagnostic use. Each human donor/unit used in the preparation of the calibrator has been tested by an FDA licensed method/test and found to be negative or non-reactive for the presence of HBsAg, Anti-HCV, NAT testing for HIV-1, HCV (RNA) and HIV-1/2. Each unit is also negative by a serological test for Syphilis (RPR or STS). Because no test method can offer complete assurance that infectious agents are absent, this material should be handled as potentially infectious. When handling or disposing of tubes follow precautions for patient specimens as specified in the OSHA Bloodborne Pathogen Rule (29 CFR Part 1910, 1030) or other equivalent biosafety procedures.
- When handling the calibrator, wear protective equipment such as disposable gloves to protect yourself from infection.
- Do not swallow the calibrator. It may be infectious.

- NOTE**
- Do not use the calibrator if it is past the expiration date on the calibrator or if the calibrator is stored under unspecified conditions.
 - Store the calibrator in a pharmaceutical refrigerator.
 - Store the calibrator at temperatures between 2 and 8°C (36 and 46°F).
 - Do not freeze the calibrator.
 - Accurate values cannot be obtained if the calibrator is used in a non-mechanical hematological testing.

Using the Calibrator

- NOTE**
- Refer to the analyzer operator's manual for details on measurement.
 - Do not use the calibrator if deterioration is suspected. After mixing, the calibrator should be similar in appearance to fresh whole blood. In unmixed tubes, the supernatant may appear cloudy and reddish; this is normal and does not indicate deterioration. Other discoloration, very dark red supernatant or unacceptable results may indicate deterioration.
 - Confirm that the lot number of the calibrator matches the lot number on the assay sheet.

- When using the MEK-9100 or MEK-9200 automated hematology analyzer, use the handy barcode reader provided with the hematology analyzer to read the 2D barcodes on the assay sheet of a Nihon Kohden genuine calibrator. The assay values are entered into the hematology analyzer.
- When using the MEK-1301 or MEK-1302 automated hematology analyzer, MEK-1303 automated hematology and clinical chemistry analyzer or MEK-1305 automated hematology and ESR analyzer, read the barcodes printed on the "Barcode Entry" area on the assay sheet of the Nihon Kohden genuine calibrator with the connected barcode reader.

Frequency of Performing Calibration

Perform calibration in the following cases.

- Calibration is required as a result of regular quality control.
- Calibration is required as a result of quality control performed after maintenance.
- Other cases where calibration is required.

Calibrator Used


MEK-CAL
HEMATOLOGY CALIBRATOR

Measurement Principle

Refer to the analyzer operator's manual.

Procedure

1. Take the calibrator out of the pharmaceutical refrigerator and allow it to warm to room temperature (15 to 30°C, 59 to 86°F) for 15 minutes before mixing.

NOTE: Do not mix the calibrator before it warms to room temperature. Doing so may cause the contents of the calibrator to ooze out the top of the cap.


2. Mix the calibrator.

NOTE


- Do not use a mechanical mixer.
- Do not mix the calibrator roughly. If it bubbles, hemolysis is caused.

- 1) Hold the tube containing the calibrator between the palms of the hands and roll it back and forth. Occasionally invert the tube so that the calibrator is mixed thoroughly.

Roll the tube back and forth

Invert the tube

- 2) Continue to mix in the manner described in step 1) until the red cells are completely suspended. Tubes stored for a long time may require extra mixing.
3. Gently invert the tube until there is no clot of blood on the bottom before measurement.

NOTE: If mixing is incomplete, measurement results may be inaccurate.

4. Measure the calibrator as instructed in the analyzer operator's manual.

NOTE: If the tube was opened for the measurement, clean residual material from the cap and tube rim with a lint-free tissue after the measurement, then close the cap tightly.

5. Return the calibrator to the pharmaceutical refrigerator.

NOTE: Leaving the calibrator at room temperature for an extended period of time may affect the results of subsequent measurements.

Technical Information

Metrological Traceability of Values Assigned to the Calibrator

For each parameter, the internationally recognized reference method is:

- 1) WBC and RBC: ICSH 1988: The assignment of values to fresh blood used for calibrating automated blood cell counters. Clin Lab Hematol.1988; 10:203-212
- 2) HGB: CLSI H15-A3
- 3) HCT: CLSI H7-A3
- 4) PLT: International Council for Standardization in Hematology Expert Panel on Cytometry and International Society of Laboratory Hematology Task Force on Platelet Counting. Platelet counting by RBC/platelet ratio method. A reference method. Am Journal of Clinical Pathology 115:460-464 2001

Assay Values

- Refer to the assay values for the specific instrument model on the provided assay sheet.
- Assay values are determined on well-maintained, properly calibrated instruments using the conformed reagents. Reagent differences, maintenance, and operation technique may contribute to inter-laboratory variation.
- Assigned values are presented as a Mean. The Mean is derived from replicate testing on instruments operated and maintained according to the manufacturer's instructions.

NOTE: RDW-CV is displayed as RDW on the MEK-6400, MEK-6410, MEK-6420, MEK-6450, MEK-7222 and MEK-8222 automatic hematology analyzers.

Limitations

The performance of the calibrator is assured only if it is properly stored and used as described in this operator's manual. Incomplete mixing of the tube prior to use invalidates both the sample withdrawn and any remaining material in the tube.

Composition

The calibrator contains human erythrocytes, simulated leukocytes and mammalian platelets suspended in a plasma-like fluid with preservatives.

Sterilization Method

The calibrator is not intended to be sterilized or kept in a sterile environment.

Mathematical Approach upon which the Calculation of the Analytical Result is Made

The following parameters are calculated as follows ¹:

$$MCV = HCT / RBC$$

$$MCH = HGB / RBC$$

$$MCHC = HGB / HCT$$

$$MPV = PCT / PLT^2$$

¹ Without unit correction

² PCT: Platelet Crit

NOTE: RBC indices are calculated by rounding off to the whole number or the nearest tenth.

Environmental Conditions

Storage and Transport Environment

Temperature: 2 to 8°C (36 to 46°F)

Usage Environment

Temperature: 15 to 30°C (59 to 86°F)

Expiration Date

The expiration date is shown on the tube label and package.

Shelf Life Date After Opening

7 days (within the expiration date)

NOTE: Record the opening date at first use.

Package and Catalog Number

Model	Qty	Catalog Number
MEK-CAL	MEK-CAL (2.0 to 2.4 mL) × 1	CAL01
	MEK-CAL (2.0 to 2.4 mL) × 2	CAL02
	MEK-CAL (2.0 to 2.4 mL) × 3	CAL03
	MEK-CAL (2.0 to 2.4 mL) × 4	CAL04
	MEK-CAL (2.0 to 2.4 mL) × 5	CAL05
	MEK-CAL (2.0 to 2.4 mL) × 6	CAL06
	MEK-CAL (2.0 to 2.4 mL) × 7	CAL07
	MEK-CAL (2.0 to 2.4 mL) × 8	CAL08
	MEK-CAL (2.0 to 2.4 mL) × 9	CAL09
	MEK-CAL (2.0 to 2.4 mL) × 10	CAL10
	MEK-CAL (2.0 to 2.4 mL) × 11	CAL11
	MEK-CAL (2.0 to 2.4 mL) × 12	CAL12

Disposal

WARNING

Dispose of the calibrator according to your local laws and your facility's guidelines (including incineration, melt treatment, sterilization, disinfection and request for waste disposal) for disposing of infectious medical waste. Otherwise, it may affect the environment. If there is a possibility that the calibrator may have been contaminated with infection, it may cause infection.

When disposing of the calibrator, such as when the expiration date is past, follow the instructions on the SDS of the calibrator.

Revision History

Edition	Date	Details	Code Number
1st Edition	13 Jun 2007	Initial issue	0604-021476
7th Edition	19 Jan 2024	Corrections	0614-908008A
8th Edition	22 Oct 2024	Corrections	0614-908008B

- NOTE
- The code number of this manual was changed from 0604-021476D to 0614-908008 when the manual was updated from 5th Edition to 6th Edition.
 - Changes made in the most recent edition are indicated by a bar in the left margin of each page.

Note for users in the territory of the EEA and Switzerland:
Any serious incident that has occurred in relation to the device must be reported to the European Representative designated by the manufacturer and the Competent Authority of the Member State of the EEA and Switzerland in which the user and/or patient is established.

Copyright Notice

The entire contents of this manual are copyrighted by Nihon Kohden. All rights are reserved.

 Manufacturer
NIHON KOHDEN CORPORATION
1-31-4 Nishiochiai, Shinjuku-ku,
Tokyo 161-8560, Japan
Phone +81 3-5996-8041
<https://www.nihonkohden.com/>

 European Representative
NIHON KOHDEN EUROPE GmbH
Raiffeisenstrasse 10, 61191 Rosbach, Germany
Phone +49 6003-827-0 Fax +49 6003-827-599
UK Responsible Person
NIHON KOHDEN UK LTD.
Unit 3, Heyworth Business Park,
Old Portsmouth Road, Peasmarsh,
Guildford, Surrey, GU3 1AF, UK
Phone +44 14-8333-1328

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Za svaki parametar, međunarodno priznata referentna metoda je:

- 1) WBC i RBC: ICSH 1988: The assignment of values to fresh blood used for calibrating automated blood cell counters. Clin Lab Hematol. 1988.; 10: 203-212
- 2) HGB: CLSI H15-A3
- 3) HCT: CLSI H7-A3
- 4) PLT: International Council for Standardization in Hematology Expert Panel on Cytometry and International Society of Laboratory Hematology Task Force on Platelet Counting. Platelet counting by RBC/platelet ratio method. A reference method. Am Journal of Clinical Pathology 115:460-464 2001

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- 1) WBC a RBC: ICSH 1988: The assignment of values to fresh blood used for calibrating automated blood cell counters. Clin Lab Hematol.1988; 10:203-212
- 2) HGB: CLSI H15-A3
- 3) HCT: CLSI H7-A3
- 4) PLT: International Council for Standardization in Hematology Expert Panel on Cytometry and International Society of Laboratory Hematology Task Force on Platelet Counting. Platelet counting by RBC/platelet ratio method. A reference method. Am Journal of Clinical Pathology 115:460-464 2001

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- 1) WBC ja RBC: ICSH 1988: The assignment of values to fresh blood used for calibrating automated blood cell counters. Clin Lab Hematol.1988; 10:203-212
- 2) HGB: CLSI H15-A3
- 3) HCT: CLSI H7-A3
- 4) PLT: International Council for Standardization in Hematology Expert Panel on Cytometry and International Society of Laboratory Hematology Task Force on Platelet Counting. Platelet counting by RBC/platelet ratio method. A reference method. Am Journal of Clinical Pathology 115:460-464 2001

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- 1) WBC et RBC: ICSH 1988: The assignment of values to fresh blood used for calibrating automated blood cell counters. Clin Lab Hematol.1988; 10:203-212
- 2) HGB : CLSI H15-A3
- 3) HCT : CLSI H7-A3
- 4) PLT : International Council for Standardization in Hematology Expert Panel on Cytometry and International Society of Laboratory Hematology Task Force on Platelet Counting. Platelet counting by RBC/platelet ratio method. A reference method. Am Journal of Clinical Pathology 115:460-464 2001

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- 1) WBC und RBC: ICSH 1988: The assignment of values to fresh blood used for calibrating automated blood cell counters. Clin Lab Hematol.1988; 10:203-212
- 2) HGB: CLSI H15-A3
- 3) HCT: CLSI H7-A3
- 4) PLT: International Council for Standardization in Hematology Expert Panel on Cytometry and International Society of Laboratory Hematology Task Force on Platelet Counting. Platelet counting by RBC/platelet ratio method. A reference method. Am Journal of Clinical Pathology 115:460-464 2001

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- 1) WBC και RBC: ICSH 1988: The assignment of values to fresh blood used for calibrating automated blood cell counters. Clin Lab Hematol.1988; 10:203-212
- 2) HGB: CLSI H15-A3
- 3) HCT: CLSI H7-A3
- 4) PLT: International Council for Standardization in Hematology Expert Panel on Cytometry and International Society of Laboratory Hematology Task Force on Platelet Counting. Platelet counting by RBC/platelet ratio method. A reference method. Am Journal of Clinical Pathology 115:460-464 2001

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- 1) WBC e RBC: ICSH 1988: The assignment of values to fresh blood used for calibrating automated blood cell counters. Clin Lab Hematol.1988; 10:203-212
- 2) HGB: CLSI H15-A3
- 3) HCT: CLSI H7-A3
- 4) PLT: International Council for Standardization in Hematology Expert Panel on Cytometry and International Society of Laboratory Hematology Task Force on Platelet Counting. Platelet counting by RBC/platelet ratio method. A reference method. Am Journal of Clinical Pathology 115:460-464 2001

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- 1) WBC и RBC: ICSH 1988: The assignment of values to fresh blood used for calibrating automated blood cell counters. Clin Lab Hematol.1988; 10:203-212
- 2) HGB: CLSI H15-A3
- 3) HCT: CLSI H7-A3
- 4) PLT: International Council for Standardization in Hematology Expert Panel on Cytometry and International Society of Laboratory Hematology Task Force on Platelet Counting. Platelet counting by RBC/platelet ratio method. A reference method. Am Journal of Clinical Pathology 115:460-464 2001

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- 1) WBC in RBC: ICSH 1988: The assignment of values to fresh blood used for calibrating automated blood cell counters. Clin Lab Hematol.1988; 10:203-212
- 2) HGB: CLSI H15-A3
- 3) HCT: CLSI H7-A3

4) PLT:

International Council for Standardization in Hematology Expert Panel on Cytometry and International Society of Laboratory Hematology Task Force on Platelet Counting. Platelet counting by RBC/platelet ratio method. A reference method. Am Journal of Clinical Pathology 115:460-464 2001

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- 1) WBC y RBC: ICSH 1988: The assignment of values to fresh blood used for calibrating automated blood cell counters. Clin Lab Hematol.1988; 10:203-212
- 2) HGB: CLSI H15-A3
- 3) HCT: CLSI H7-A3
- 4) PLT: International Council for Standardization in Hematology Expert Panel on Cytometry and International Society of Laboratory Hematology Task Force on Platelet Counting. Platelet counting by RBC/platelet ratio method. A reference method. Am Journal of Clinical Pathology 115:460-464 2001

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