


MK-3CH HEMATOLOGY CONTROL+CRP 2.5mL

ASSAY SHEET

LOT	MK3C248H
	2024-11-05
OPEN	14DAYS

MEK-1303 Barcode Entry

If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.

1



2



3



		MEK-1303	RANGE
WBC	10 ³ /μL	20.35	± 2.20
RBC	10 ⁶ /μL	5.87	± 0.25
HGB	g/dL	19.32	± 0.70
HCT	%	62.9	± 4.0
MCV	fL	107.1	± 7.0
MCH	pg	32.9	± 2.8
MCHC	g/dL	30.7	± 3.0
PLT	10 ³ /μL	541.5	± 70.0
MPV	fL	8.6	± 3.0
RDW-CV	%	16.5	± 3.0
LY%	%	13.76	± 10.00
MO%	%	2.57	± 5.00
GR%	%	83.67	± 10.00
LY	10 ³ /μL	2.80	± 1.80
MO	10 ³ /μL	0.52	± 0.70
GR	10 ³ /μL	17.03	± 2.00
CRP	Precision	≤ 10.0%	
		(52.90-98.30 mg/L)	

Manual Registration	
Lot	248
Level	High
Control	MK-3C


<https://www.nihonkohden.com/spi/index.html>

NIHON KOHDEN CORPORATION

 1-31-4 NISHIOCHIAI, SHINJUKU-KU
 TOKYO 161-8560 JAPAN

0694-901767G