


**MK-RE1 HEMATOLOGY CONTROL**  
for reticulocyte **3.0mL**

**ASSAY SHEET**

LOT	MKRE249-1
	2024-11-05
OPEN	14 DAYS

		MEK-9200	RANGE
RET%	%	1.57	±1.43
RET	10 <sup>6</sup> /μL	0.0550	±0.0501

**2D Barcode Entry**



• If RBC value is under 3.00 [ $\times 10^6/\mu\text{L}$ ], please refrain from using.  
It may have deteriorated or hemolyzed, which affects the measurement accuracy.



<https://www.nihonkohden.com/spi/index.html>

Manual Registration						
Lot	249					
Level	1					
Expiration	11	M	5	D	2024	Y
Control	MKRE					



**NIHON KOHDEN CORPORATION**  
1-31-4 NISHIOCHIAI, SHINJUKU-KU  
TOKYO 161-8560 JAPAN