


MEK-CAL HEMATOLOGY CALIBRATOR 2.0mL

ASSAY SHEET

LOT	PLUS24X
	2024-11-05
OPEN	7DAYS

MEK-9100,9200
2D Barcode Entry



MEK-1301,1302,1303,1305
Barcode Entry




		MEK -6400 -6410 -6420 -6450	MEK -6500 -6510 -6550	MEK -7222	MEK -7300	MEK -8222	MEK -9100 -9200	MEK -1301 -1302 -1303 -1305
WBC	10 ³ /μL	9.1	9.1	9.3	9.2	9.0	9.39	8.99
RBC	10 ⁶ /μL	4.62	4.62	4.59	4.62	4.58	4.57	4.60
HGB	g/dL	13.0	13.1	13.0	12.9	12.8	13.12	13.11
HCT	%	38.9	38.9	40.7	39.5	38.4	43.9	44.2
MCV	fL	84.2	84.2	88.6	85.4	83.8	96.0	96.0
MCH	pg	28.1	28.4	28.3	27.9	27.9	28.7	28.5
MCHC	g/dL	33.4	33.7	31.9	32.7	33.3	29.9	29.7
PLT	10 ³ /μL	238	238	231	234	233	250.4	260.8
MPV	fL	6.9	8.9	6.9	7.7	6.9	9.4	9.4
RDW-CV*	%	16.2	16.2	14.9	16.3	15.2	16.4	16.3

MEK	-7222
	-7300
	-8222
	-9100
	-9200
FS	163
FL	89
SD	106
FS THR.	100

* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450,7222,8222.

 <https://www.nihonkohden.com/spi/index.html>

 **NIHON KOHDEN CORPORATION**
1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN



 
0123
0694-901780T