



MEK-3DH HEMATOLOGY CONTROL 2.0mL

ASSAY SHEET

LOT	B24YH
	2025-02-05
OPEN	14DAYS

MEK-1301,1302,1305 Barcode Entry

The barcodes may be read in any order. If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.



		RANGE	MEK -6400 -6410 -6420 -6450	MEK -6500 -6510 -6550	MEK -1301 -1302 -1305
WBC	10 ³ /μL	±2.2	20.6	20.6	19.89
RBC	10 ⁶ /μL	±0.25	5.97	5.97	5.97
HGB	g/dL	±0.7	19.2	19.2	19.23
HCT	%	±4.0	55.0	55.0	62.4
MCV	fL	±9.0	92.1	92.1	104.6
MCH	pg	±2.8	32.2	32.2	32.2
MCHC	g/dL	±3.0	34.9	34.9	30.8
PLT	10 ³ /μL	±70	558	558	559.5
MPV	fL	±3.0	5.6	8.1	8.5
RDW-CV *	%	±3.0	15.9	15.9	17.7
LY%	%	±10.0	21.4	21.4	17.49
MO%	%	±5.0	1.4	1.4	1.14
GR%	%	±10.0	77.2	77.2	81.37
LY	10 ³ /μL	±1.8	4.4	4.4	3.48
MO	10 ³ /μL	±0.7	0.3	0.3	0.23
GR	10 ³ /μL	±2.0	15.9	15.9	16.18

Manual Registration	
Lot	24Y
Level	High
Control	MEK-3D

* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450.



<https://www.nihonkohden.com/spi/index.html>

NIHON KOHDEN CORPORATION
 1-31-4 NISHIOCHIAI, SHINJUKU-KU
 TOKYO 161-8560 JAPAN

