


<b>MEK-3DL</b>	HEMATOLOGY CONTROL	2.0mL
----------------	--------------------	-------

**ASSAY SHEET**

LOT	B24YL
	2025-02-05
OPEN	14DAYS

<b>MEK-1301,1302,1305 Barcode Entry</b>
---

The barcodes may be read in any order.  
 If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.

1



2



3



		RANGE	MEK -6400 -6410 -6420 -6450	MEK -6500 -6510 -6550	MEK -1301 -1302 -1305
<b>WBC</b>	10 <sup>3</sup> /μL	±0.6	2.1	2.1	2.03
<b>RBC</b>	10 <sup>6</sup> /μL	±0.15	2.47	2.47	2.40
<b>HGB</b>	g/dL	±0.4	6.0	6.1	6.05
<b>HCT</b>	%	±2.0	18.4	18.4	20.4
<b>MCV</b>	fL	±9.0	74.5	74.5	84.8
<b>MCH</b>	pg	±2.4	24.3	24.7	25.2
<b>MCHC</b>	g/dL	±3.0	32.6	33.2	29.7
<b>PLT</b>	10 <sup>3</sup> /μL	±30	61	61	69.3
<b>MPV</b>	fL	±3.0	5.9	8.3	8.2
<b>RDW-CV *</b>	%	±3.0	17.1	17.1	15.7
<b>LY%</b>	%	±10.0	62.4	62.4	62.71
<b>MO%</b>	%	±5.0	6.2	6.2	3.54
<b>GR%</b>	%	±10.0	31.4	31.4	33.75
<b>LY</b>	10 <sup>3</sup> /μL	±0.5	1.3	1.3	1.27
<b>MO</b>	10 <sup>3</sup> /μL	±0.2	0.1	0.1	0.07
<b>GR</b>	10 <sup>3</sup> /μL	±0.5	0.7	0.7	0.69

<b>MEK-1305</b>
QC ESR
44±17 mm

Manual Registration	
Lot	<b>24Y</b>
Level	<b>Low</b>
Control	<b>MEK-3D</b>

\* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450.



<https://www.nihonkohden.com/spi/index.html>



**NIHON KOHDEN CORPORATION**  
 1-31-4 NISHIOCHIAI, SHINJUKU-KU  
 TOKYO 161-8560 JAPAN



0694-901684R