


<b>MEK-3DN</b>	<b>HEMATOLOGY CONTROL</b>	<b>2.0mL</b>
----------------	---------------------------	--------------

**ASSAY SHEET**

<b>LOT</b>	B24YN
	2025-02-05
<b>OPEN</b>	14DAYS

<b>MEK-1301,1302,1305 Barcode Entry</b>
---

The barcodes may be read in any order.  
 If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.

1



2



3



		RANGE	MEK -6400	MEK -6500	MEK -1301
			-6410	-6510	-1302
			-6420	-6550	-1305
			-6450		
<b>WBC</b>	10 <sup>3</sup> /μL	±0.8	7.6	7.6	7.58
<b>RBC</b>	10 <sup>6</sup> /μL	±0.20	4.83	4.83	4.82
<b>HGB</b>	g/dL	±0.5	13.7	13.9	13.92
<b>HCT</b>	%	±3.0	40.9	40.9	46.1
<b>MCV</b>	fL	±9.0	84.7	84.7	95.6
<b>MCH</b>	pg	±2.4	28.4	28.8	28.9
<b>MCHC</b>	g/dL	±3.0	33.5	34.0	30.2
<b>PLT</b>	10 <sup>3</sup> /μL	±50	257	257	278.7
<b>MPV</b>	fL	±2.0	5.6	8.1	8.5
<b>RDW-CV *</b>	%	±3.0	16.5	16.5	16.9
<b>LY%</b>	%	±10.0	38.3	38.3	34.94
<b>MO%</b>	%	±5.0	4.0	4.0	2.00
<b>GR%</b>	%	±10.0	57.7	57.7	63.06
<b>LY</b>	10 <sup>3</sup> /μL	±1.0	2.9	2.9	2.65
<b>MO</b>	10 <sup>3</sup> /μL	±0.4	0.3	0.3	0.15
<b>GR</b>	10 <sup>3</sup> /μL	±1.0	4.4	4.4	4.78

<b>MEK-1305</b>
QC ESR
11±5 mm

Manual Registration	
Lot	24Y
Level	Normal
Control	MEK-3D

\* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450.



<https://www.nihonkohden.com/spi/index.html>



**NIHON KOHDEN CORPORATION**

1-31-4 NISHIOCHIAI, SHINJUKU-KU  
 TOKYO 161-8560 JAPAN



MEK-3DN



B24YN



2025-02-05



0694-901685R