


**MEK-5DH** HEMATOLOGY CONTROL  
for WBC 5 part differential 3.0mL

**ASSAY SHEET**

LOT	MEK5D24YH
	2025-02-05
OPEN	14DAYS

MEK-9100,9200  
2D Barcode Entry



		RANGE	MEK -7222	MEK -7300	MEK -8222	MEK -9100 -9200
WBC	10 <sup>3</sup> /μL	±2.2	25.3	25.2	24.7	26.09
NE%	%	±10.0	58.4	58.4	58.4	61.98
LY%	%	±10.0	17.9	17.9	17.9	19.37
MO%	%	±5.0	17.0	17.0	17.0	10.73
EO%	%	±5.0	5.6	5.6	5.6	7.06
BA%	%	±4.0	1.1	1.1	1.1	0.86
NE	10 <sup>3</sup> /μL	±2.8	14.8	14.7	14.4	16.18
LY	10 <sup>3</sup> /μL	±2.8	4.5	4.5	4.4	5.05
MO	10 <sup>3</sup> /μL	±1.0	4.3	4.3	4.2	2.80
EO	10 <sup>3</sup> /μL	±1.0	1.4	1.4	1.4	1.84
BA	10 <sup>3</sup> /μL	±0.8	0.3	0.3	0.3	0.22
RBC	10 <sup>6</sup> /μL	±0.25	5.38	5.34	5.33	5.39
HGB	g/dL	±0.7	16.6	16.4	16.5	16.51
HCT	%	±4.0	52.5	50.9	49.5	58.3
MCV	fL	±8.0	97.5	95.4	92.9	108.2
MCH	pg	±2.8	30.9	30.7	31.0	30.6
MCHC	g/dL	±3.0	31.6	32.2	33.3	28.3
RDW-CV *	%	±4.0	13.8	14.7	14.0	15.3
PLT	10 <sup>3</sup> /μL	±70	563	554	562	564.0
MPV	fL	±4.0	9.3	8.9	9.5	11.0

Manual Registration	
Lot	24Y
Level	High
Control	MEK5D

\* RDW-CV is displayed as RDW on MEK-7222,8222.



<https://www.nihonkohden.com/spi/index.html>

**NIHON KOHDEN CORPORATION**  
1-31-4 NISHIOCHIAI, SHINJUKU-KU  
TOKYO 161-8560 JAPAN



0694-901692Q