


MEK-5DL HEMATOLOGY CONTROL
for WBC 5 part differential 3.0mL

ASSAY SHEET

LOT	MEK5D24YL
	2025-02-05
OPEN	14DAYS

MEK-9100,9200
2D Barcode Entry



		RANGE	MEK -7222	MEK -7300	MEK -8222	MEK -9100 -9200
WBC	10 ³ /μL	±0.6	2.9	2.9	2.9	2.96
NE%	%	±15.0	31.0	31.0	31.0	30.85
LY%	%	±15.0	48.3	48.3	48.3	51.75
MO%	%	±5.0	13.9	13.9	13.9	11.27
EO%	%	±5.0	4.0	4.0	4.0	4.54
BA%	%	±4.0	2.8	2.8	2.8	1.59
NE	10 ³ /μL	±0.5	0.9	0.9	0.9	0.92
LY	10 ³ /μL	±0.5	1.4	1.4	1.4	1.53
MO	10 ³ /μL	±0.2	0.4	0.4	0.4	0.33
EO	10 ³ /μL	±0.2	0.1	0.1	0.1	0.13
BA	10 ³ /μL	±0.2	0.1	0.1	0.1	0.05
RBC	10 ⁶ /μL	±0.15	2.63	2.61	2.63	2.51
HGB	g/dL	±0.4	6.3	6.4	6.2	6.44
HCT	%	±2.0	21.6	20.8	20.4	22.3
MCV	fL	±7.0	82.2	79.7	77.6	89.0
MCH	pg	±2.4	24.0	24.5	23.6	25.7
MCHC	g/dL	±3.0	29.2	30.8	30.4	28.9
RDW-CV *	%	±4.0	14.8	16.1	15.1	14.9
PLT	10 ³ /μL	±30	51	57	51	51.9
MPV	fL	±3.0	4.1	5.2	4.4	6.9

Manual Registration	
Lot	24Y
Level	Low
Control	MEK5D

* RDW-CV is displayed as RDW on MEK-7222,8222.



<https://www.nihonkohden.com/spi/index.html>

NIHON KOHDEN CORPORATION
1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN



0694-901690Q