


MEK-5DN HEMATOLOGY CONTROL
for WBC 5 part differential 3.0mL

ASSAY SHEET

LOT	MEK5D24YN
	2025-02-05
OPEN	14DAYS

MEK-9100,9200
2D Barcode Entry



		RANGE	MEK -7222	MEK -7300	MEK -8222	MEK -9100 -9200
WBC	10 ³ /μL	±0.8	8.3	8.2	8.1	8.34
NE%	%	±15.0	43.8	43.8	43.8	43.25
LY%	%	±15.0	37.8	37.8	37.8	42.52
MO%	%	±5.0	13.6	13.6	13.6	10.20
EO%	%	±5.0	1.4	1.4	1.4	1.43
BA%	%	±4.0	3.4	3.4	3.4	2.60
NE	10 ³ /μL	±1.5	3.7	3.6	3.5	3.60
LY	10 ³ /μL	±1.5	3.1	3.1	3.1	3.55
MO	10 ³ /μL	±0.5	1.1	1.1	1.1	0.85
EO	10 ³ /μL	±0.5	0.1	0.1	0.1	0.12
BA	10 ³ /μL	±0.4	0.3	0.3	0.3	0.22
RBC	10 ⁶ /μL	±0.20	4.59	4.52	4.57	4.50
HGB	g/dL	±0.5	13.2	13.2	13.0	13.35
HCT	%	±3.0	42.0	40.3	39.6	45.4
MCV	fL	±6.0	91.6	89.1	86.7	100.8
MCH	pg	±2.8	28.8	29.2	28.4	29.7
MCHC	g/dL	±3.0	31.4	32.8	32.8	29.4
RDW-CV *	%	±4.0	14.3	15.3	14.5	15.1
PLT	10 ³ /μL	±50	269	286	263	260.8
MPV	fL	±2.0	3.3	5.1	3.5	6.8

Manual Registration	
Lot	24Y
Level	Normal
Control	MEK5D

* RDW-CV is displayed as RDW on MEK-7222,8222.



<https://www.nihonkohden.com/spi/index.html>

NIHON KOHDEN CORPORATION
1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN



0694-901691Q