


MK-3CH HEMATOLOGY CONTROL+CRP 2.5mL

ASSAY SHEET

LOT	MK3C24YH
	2025-02-05
OPEN	14DAYS

MEK-1303 Barcode Entry

If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.

1



2



3



		MEK-1303	RANGE
WBC	10 ³ /μL	19.80	± 2.20
RBC	10 ⁶ /μL	5.95	± 0.25
HGB	g/dL	19.17	± 0.70
HCT	%	62.7	± 4.0
MCV	fL	105.4	± 9.0
MCH	pg	32.2	± 2.8
MCHC	g/dL	30.6	± 3.0
PLT	10 ³ /μL	556.1	± 70.0
MPV	fL	8.5	± 3.0
RDW-CV	%	17.1	± 3.0
LY%	%	17.74	± 10.00
MO%	%	1.29	± 5.00
GR%	%	80.97	± 10.00
LY	10 ³ /μL	3.51	± 1.80
MO	10 ³ /μL	0.26	± 0.70
GR	10 ³ /μL	16.03	± 2.00
CRP	Precision	≤ 10.0%	
		(50.20-93.20 mg/L)	

Manual Registration	
Lot	24Y
Level	High
Control	MK-3C


<https://www.nihonkohden.com/spi/index.html>

 **NIHON KOHDEN CORPORATION**
 1-31-4 NISHIOCHIAI, SHINJUKU-KU
 TOKYO 161-8560 JAPAN



0694-901767H