

MEK-CAL HEMATOLOGY CALIBRATOR 2.0mL

ASSAY SHEET

LOT	PLUS24Z
	2025-01-05
OPEN	7DAYS

MEK-9100,9200
2D Barcode Entry



MEK-1301,1302,1303,1305
Barcode Entry



		MEK -6400 -6410 -6420 -6450	MEK -6500 -6510 -6550	MEK -7222	MEK -7300	MEK -8222	MEK -9100 -9200	MEK -1301 -1302 -1303 -1305
WBC	10 ³ /μL	8.9	8.9	9.2	9.1	8.9	9.27	8.94
RBC	10 ⁶ /μL	4.72	4.72	4.68	4.70	4.67	4.67	4.70
HGB	g/dL	13.5	13.7	13.4	13.4	13.3	13.62	13.61
HCT	%	40.1	40.1	41.7	40.4	39.3	45.3	45.6
MCV	fL	84.9	84.9	89.0	85.9	84.2	97.0	97.1
MCH	pg	28.6	29.0	28.6	28.5	28.5	29.2	29.0
MCHC	g/dL	33.7	34.2	32.1	33.2	33.8	30.1	29.8
PLT	10 ³ /μL	241	241	238	241	241	255.8	261.7
MPV	fL	6.0	8.4	6.0	7.2	6.0	8.6	8.7
RDW-CV*	%	16.4	16.4	15.1	16.6	15.5	16.8	16.7

MEK	-7222
	-7300
	-8222
	-9100
	-9200
FS	172
FL	86
SD	114
FS THR.	100

* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450,7222,8222.



<https://www.nihonkohden.com/spi/index.html>

NIHON KOHDEN CORPORATION
1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN



0694-901780V