



MEK-CAL HEMATOLOGY CALIBRATOR 2.0mL

ASSAY SHEET

| | |
|---|------------|
| LOT | PLUS252 |
|  | 2025-03-05 |
| OPEN | 7DAYS |

MEK-9100,9200
2D Barcode Entry



MEK-1301,1302,1303,1305
Barcode Entry



| | | MEK -6400 -6410 -6420 -6450 | MEK -6500 -6510 -6550 | MEK -7222 | MEK -7300 | MEK -8222 | MEK -9100 -9200 | MEK -1301 -1302 -1303 -1305 |
|----------|---------------------|---|--------------------------------|--------------|--------------|--------------|-----------------------|---|
| WBC | 10 ³ /μL | 9.3 | 9.3 | 9.6 | 9.5 | 9.3 | 9.69 | 9.30 |
| RBC | 10 ⁶ /μL | 4.70 | 4.70 | 4.68 | 4.70 | 4.67 | 4.68 | 4.70 |
| HGB | g/dL | 13.7 | 13.8 | 13.6 | 13.5 | 13.4 | 13.75 | 13.66 |
| HCT | % | 40.0 | 40.0 | 41.9 | 40.7 | 39.6 | 46.1 | 46.4 |
| MCV | fL | 85.2 | 85.2 | 89.5 | 86.5 | 84.7 | 98.4 | 98.7 |
| MCH | pg | 29.1 | 29.4 | 29.1 | 28.7 | 28.7 | 29.4 | 29.1 |
| MCHC | g/dL | 34.3 | 34.5 | 32.5 | 33.2 | 33.8 | 29.8 | 29.4 |
| PLT | 10 ³ /μL | 242 | 242 | 234 | 237 | 236 | 253.2 | 263.0 |
| MPV | fL | 6.4 | 8.6 | 6.4 | 7.4 | 6.4 | 9.0 | 9.0 |
| RDW-CV * | % | 15.5 | 15.5 | 14.5 | 15.6 | 14.7 | 15.6 | 15.5 |

| | |
|---------|-------|
| MEK | -7222 |
| | -7300 |
| | -8222 |
| | -9100 |
| | -9200 |
| FS | 168 |
| FL | 87 |
| SD | 105 |
| FS THR. | 100 |

* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450,7222,8222.



<https://www.nihonkohden.com/spi/index.html>

NIHON KOHDEN CORPORATION
1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN

