


MEK-CAL HEMATOLOGY CALIBRATOR 2.0mL

ASSAY SHEET

LOT	PLUS253
	2025-04-05
OPEN	7DAYS

MEK-9100,9200
2D Barcode Entry



MEK-1301,1302,1303,1305
Barcode Entry



		MEK -6400 -6410 -6420 -6450	MEK -6500 -6510 -6550	MEK -7222	MEK -7300	MEK -8222	MEK -9100 -9200	MEK -1301 -1302 -1303 -1305
WBC	10 ³ /μL	9.0	9.0	9.2	9.1	9.0	9.26	8.88
RBC	10 ⁶ /μL	4.59	4.59	4.56	4.60	4.55	4.53	4.61
HGB	g/dL	13.4	13.6	13.3	13.3	13.2	13.56	13.47
HCT	%	39.7	39.7	41.6	40.4	39.3	44.9	45.9
MCV	fL	86.5	86.5	91.3	87.8	86.4	99.2	99.5
MCH	pg	29.2	29.6	29.2	28.9	29.0	29.9	29.2
MCHC	g/dL	33.8	34.3	32.0	32.9	33.6	30.2	29.3
PLT	10 ³ /μL	232	232	224	227	226	242.8	253.1
MPV	fL	6.1	8.5	6.1	7.3	6.1	8.9	8.9
RDW-CV*	%	15.6	15.6	14.5	15.7	14.8	15.6	15.5

MEK	-7222
	-7300
	-8222
	-9100
	-9200
FS	167
FL	88
SD	103
FS THR.	100

* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450,7222,8222.



<https://www.nihonkohden.com/spi/index.html>



NIHON KOHDEN CORPORATION
1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN



MEK-CAL



PLUS253



2025-04-05



0694-901780Y