




MK-RE2

HEMATOLOGY CONTROL
for reticulocyte


3.0mL

ASSAY SHEET

 LOT	MKRE253-2
	2025-05-05
 OPEN	14 DAYS

		MEK-9200	RANGE
RET%	%	4.46	±2.54
RET	10⁶/μL	0.1553	±0.0885

2D Barcode Entry



- If RBC value is under 3.00 [$\times 10^6/\mu\text{L}$], please refrain from using.
It may have deteriorated or hemolyzed, which affects the measurement accuracy.



<https://www.nihonkohden.com/spi/index.html>

Manual Registration

Lot	253					
Level	2					
Expiration	5	M	5	D	2025	Y
Control	MKRE					



NIHON KOHDEN CORPORATION
1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN