


MEK-CAL HEMATOLOGY CALIBRATOR 2.0mL

ASSAY SHEET

LOT	PLUS254
	2025-05-05
OPEN	7DAYS

		MEK -6400	MEK -6500	MEK -7222	MEK -7300	MEK -8222	MEK -9100	MEK -1301
		-6410	-6510				-9200	-1302
		-6420	-6550					-1303
		-6450						-1305
WBC	10 ³ /μL	9.0	9.0	9.2	9.1	8.9	9.35	8.97
RBC	10 ⁶ /μL	4.67	4.67	4.66	4.68	4.65	4.65	4.66
HGB	g/dL	13.4	13.6	13.3	13.3	13.2	13.56	13.47
HCT	%	39.9	39.9	41.8	40.5	39.4	45.5	45.5
MCV	fL	85.4	85.4	89.7	86.6	84.8	97.9	97.7
MCH	pg	28.7	29.1	28.5	28.4	28.4	29.2	28.9
MCHC	g/dL	33.6	34.1	31.8	32.8	33.5	29.8	29.6
PLT	10 ³ /μL	223	223	220	222	222	238.5	242.9
MPV	fL	5.6	8.2	5.6	7.0	5.6	8.4	8.5
RDW-CV*	%	16.0	16.0	14.9	16.2	15.2	16.1	16.0

MEK	-7222
	-7300
	-8222
	-9100
	-9200
FS	169
FL	88
SD	108
FS THR.	100

MEK-9100,9200
2D Barcode Entry



MEK-1301,1302,1303,1305
Barcode Entry



* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450,7222,8222.



<https://www.nihonkohden.com/spi/index.html>

NIHON KOHDEN CORPORATION
1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN



MEK-CAL



PLUS254



2025-05-05



0694-901780Z