

MEK-3DH

HEMATOLOGY CONTROL

2.0mL

ASSAY SHEET

LOT	B255H
	2025-08-05
OPEN	14DAYS

MEK-1301,1302,1305 Barcode Entry

The barcodes may be read in any order.

If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.

1



2



3



		RANGE	MEK -6400 -6410 -6420 -6450	MEK -6500 -6510 -6550	MEK -1301 -1302 -1305
WBC	10 ³ /μL	±2.2	21.2	21.2	20.54
RBC	10 ⁶ /μL	±0.25	5.99	5.99	6.08
HGB	g/dL	±0.7	19.5	19.5	19.57
HCT	%	±4.0	55.6	55.6	64.6
MCV	fL	±7.0	92.9	92.9	106.3
MCH	pg	±2.8	32.6	32.6	32.2
MCHC	g/dL	±3.0	35.1	35.1	30.3
PLT	10 ³ /μL	±70	531	531	549.3
MPV	fL	±3.0	5.1	7.9	8.3
RDW-CV *	%	±3.0	16.1	16.1	17.9
LY%	%	±10.0	21.2	21.2	14.83
MO%	%	±5.0	1.1	1.1	0.55
GR%	%	±10.0	77.7	77.7	84.62
LY	10 ³ /μL	±1.8	4.5	4.5	3.05
MO	10 ³ /μL	±0.7	0.2	0.2	0.11
GR	10 ³ /μL	±2.0	16.5	16.5	17.38

* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450.

Manual Registration

Lot	255
Level	High
Control	MEK-3D


<https://www.nihonkohden.com/spi/index.html>

NIHON KOHDEN CORPORATION

1-31-4 NISHIOCHIAI, SHINJUKU-KU

TOKYO 161-8560 JAPAN



MEK-3DH



B255H



2025-08-05



0694-901686T