


| | |
|---|-------------------|
| LOT | B255N |
|  | 2025-08-05 |
| OPEN | 14DAYS |

MEK-1301,1302,1305 Barcode Entry

The barcodes may be read in any order.

If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.

1



2



3



| | | RANGE | MEK | MEK | MEK |
|-----------------|---------------------|--------------|--------------|--------------|--------------|
| | | | -6400 | -6500 | -1301 |
| | | | -6410 | -6510 | -1302 |
| | | | -6420 | -6550 | -1305 |
| | | | -6450 | | |
| WBC | 10 ³ /μL | ±0.8 | 8.0 | 8.0 | 7.82 |
| RBC | 10 ⁶ /μL | ±0.20 | 4.89 | 4.89 | 4.91 |
| HGB | g/dL | ±0.5 | 14.1 | 14.2 | 14.11 |
| HCT | % | ±3.0 | 41.6 | 41.6 | 47.5 |
| MCV | fL | ±7.0 | 85.0 | 85.0 | 96.7 |
| MCH | pg | ±2.4 | 28.8 | 29.0 | 28.7 |
| MCHC | g/dL | ±3.0 | 33.9 | 34.1 | 29.7 |
| PLT | 10 ³ /μL | ±50 | 251 | 251 | 274.6 |
| MPV | fL | ±2.0 | 5.3 | 8.0 | 8.2 |
| RDW-CV * | % | ±3.0 | 16.7 | 16.7 | 17.1 |
| LY% | % | ±10.0 | 33.0 | 33.0 | 30.49 |
| MO% | % | ±5.0 | 3.4 | 3.4 | 2.55 |
| GR% | % | ±10.0 | 63.6 | 63.6 | 66.96 |
| LY | 10 ³ /μL | ±1.0 | 2.6 | 2.6 | 2.38 |
| MO | 10 ³ /μL | ±0.4 | 0.3 | 0.3 | 0.20 |
| GR | 10 ³ /μL | ±1.0 | 5.1 | 5.1 | 5.24 |

| |
|-----------------|
| MEK-1305 |
| QC ESR |
| 8±5 mm |

| Manual Registration | |
|----------------------------|---------------|
| Lot | 255 |
| Level | Normal |
| Control | MEK-3D |

* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450.



<https://www.nihonkohden.com/spi/index.html>



NIHON KOHDEN CORPORATION
 1-31-4 NISHIOCHIAI, SHINJUKU-KU
 TOKYO 161-8560 JAPAN



MEK-3DN



B255N



2025-08-05



0123



0694-901685T