

**MK-3CH** HEMATOLOGY CONTROL+CRP

2.5mL

**ASSAY SHEET**

|                                                                                   |            |
|-----------------------------------------------------------------------------------|------------|
| LOT                                                                               | MK3C255H   |
|  | 2025-08-05 |
| OPEN                                                                              | 14DAYS     |

**MEK-1303 Barcode Entry**

If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.

1



2



3



|        |                     | MEK-1303                      | RANGE   |
|--------|---------------------|-------------------------------|---------|
| WBC    | 10 <sup>3</sup> /μL | 20.73                         | ± 2.20  |
| RBC    | 10 <sup>6</sup> /μL | 6.02                          | ± 0.25  |
| HGB    | g/dL                | 19.49                         | ± 0.70  |
| HCT    | %                   | 64.2                          | ± 4.0   |
| MCV    | fL                  | 106.6                         | ± 7.0   |
| MCH    | pg                  | 32.4                          | ± 2.8   |
| MCHC   | g/dL                | 30.4                          | ± 3.0   |
| PLT    | 10 <sup>3</sup> /μL | 543.2                         | ± 70.0  |
| MPV    | fL                  | 8.2                           | ± 3.0   |
| RDW-CV | %                   | 17.4                          | ± 3.0   |
| LY%    | %                   | 15.34                         | ± 10.00 |
| MO%    | %                   | 0.55                          | ± 5.00  |
| GR%    | %                   | 84.11                         | ± 10.00 |
| LY     | 10 <sup>3</sup> /μL | 3.18                          | ± 1.80  |
| MO     | 10 <sup>3</sup> /μL | 0.11                          | ± 0.70  |
| GR     | 10 <sup>3</sup> /μL | 17.44                         | ± 2.00  |
| CRP    | Precision           | ≤ 10.0%<br>(48.90-90.70 mg/L) |         |

**Manual Registration**

|         |       |
|---------|-------|
| Lot     | 255   |
| Level   | High  |
| Control | MK-3C |


<https://www.nihonkohden.com/spi/index.html>

**NIHON KOHDEN CORPORATION**

1-31-4 NISHIOCHIAI, SHINJUKU-KU  
TOKYO 161-8560 JAPAN


0694-901767J