

MK-3CL HEMATOLOGY CONTROL+CRP

2.5mL

ASSAY SHEET

LOT	MK3C255L
	2025-08-05
OPEN	14DAYS

MEK-1303 Barcode Entry

If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.

1



2



3



		MEK-1303	RANGE
WBC	10 ³ /μL	2.12	± 0.60
RBC	10 ⁶ /μL	2.32	± 0.15
HGB	g/dL	6.17	± 0.40
HCT	%	20.7	± 2.0
MCV	fL	89.3	± 9.0
MCH	pg	26.6	± 2.4
MCHC	g/dL	29.8	± 3.0
PLT	10 ³ /μL	79.3	± 30.0
MPV	fL	8.2	± 3.0
RDW-CV	%	16.7	± 3.0
LY%	%	56.67	± 10.00
MO%	%	4.40	± 5.00
GR%	%	38.93	± 10.00
LY	10 ³ /μL	1.20	± 0.50
MO	10 ³ /μL	0.09	± 0.20
GR	10 ³ /μL	0.83	± 0.50
CRP	Precision	≤ 10.0%	
		(1.60-9.40 mg/L)	

Manual Registration

Lot	255
Level	Low
Control	MK-3C


<https://www.nihonkohden.com/spi/index.html>

NIHON KOHDEN CORPORATION

1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN


0694-901765J