

MEK-CAL

HEMATOLOGY CALIBRATOR

2.0mL

ASSAY SHEET

LOT	PLUS256
	2025-07-05
OPEN	7DAYS

MEK-9100,9200
2D Barcode Entry

MEK-1301,1302,1303,1305
Barcode Entry


MEK-CAL


<https://www.nihonkohden.com/spi/index.html>


PLUS256



2025-07-05

MEK	-7222
	-7300
	-8222
	-9100
	-9200
FS	166
FL	87
SD	107
FS THR.	100

		MEK -6400 -6410 -6420 -6450	MEK -6500 -6510 -6550	MEK -7222	MEK -7300	MEK -8222	MEK -9100 -9200	MEK -1301 -1302 -1303 -1305
WBC	10 ³ /μL	9.2	9.2	9.4	9.3	9.2	9.65	9.14
RBC	10 ⁶ /μL	4.64	4.64	4.65	4.68	4.64	4.63	4.62
HGB	g/dL	13.6	13.7	13.4	13.4	13.3	13.65	13.51
HCT	%	40.1	40.1	42.3	41.0	40.0	45.9	45.7
MCV	fL	86.4	86.4	91.0	87.7	86.1	99.2	99.0
MCH	pg	29.3	29.5	28.8	28.6	28.7	29.5	29.2
MCHC	g/dL	33.9	34.2	31.7	32.7	33.3	29.7	29.6
PLT	10 ³ /μL	233	233	229	232	231	248.5	256.9
MPV	fL	6.0	8.4	6.0	7.2	6.0	8.8	8.8
RDW-CV *	%	16.1	16.1	14.9	16.3	15.2	16.2	16.1

* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450,7222,8222.


NIHON KOHDEN CORPORATION

1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN


0694-901832A