

MEK-CAL

HEMATOLOGY CALIBRATOR

2.0mL

ASSAY SHEET

LOT	PLUS257
	2025-08-05
OPEN	7DAYS

MEK-9100,9200
2D Barcode Entry

MEK-1301,1302,1303,1305
Barcode Entry

<https://www.nihonkohden.com/spi/index.html>


		MEK -6400 -6410 -6420 -6450	MEK -6500 -6510 -6550	MEK -7222	MEK -7300	MEK -8222	MEK -9100 -9200	MEK -1301 -1302 -1303 -1305
WBC	10 ³ /μL	9.1	9.1	9.3	9.2	9.0	9.52	9.04
RBC	10 ⁶ /μL	4.82	4.82	4.78	4.80	4.77	4.77	4.77
HGB	g/dL	14.1	14.3	14.0	14.0	13.9	14.22	14.08
HCT	%	41.9	41.9	43.6	42.3	41.2	47.9	47.8
MCV	fL	87.0	87.0	91.2	88.2	86.3	100.5	100.2
MCH	pg	29.3	29.7	29.3	29.2	29.1	29.8	29.5
MCHC	g/dL	33.7	34.1	32.1	33.1	33.7	29.7	29.5
PLT	10 ³ /μL	222	222	218	221	221	233.2	237.4
MPV	fL	5.8	8.3	5.8	7.1	5.8	8.5	8.6
RDW-CV *	%	15.3	15.3	14.4	15.5	14.6	15.3	15.3

MEK	-7222
	-7300
	-8222
	-9100
	-9200
FS	168
FL	86
SD	109
FS THR.	100

* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450,7222,8222.


NIHON KOHDEN CORPORATION
1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN


0694-901832B

