

MK-3CL HEMATOLOGY CONTROL+CRP 2.5mL

ASSAY SHEET

LOT	MK3C258L
	2025-11-05
OPEN	14DAYS

MEK-1303 Barcode Entry

If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.

1



2



3



		MEK-1303	RANGE
WBC	$10^3/\mu\text{L}$	2.00	± 0.60
RBC	$10^6/\mu\text{L}$	2.25	± 0.15
HGB	g/dL	6.08	± 0.40
HCT	%	20.2	± 2.0
MCV	fL	89.6	± 9.0
MCH	pg	27.0	± 2.4
MCHC	g/dL	30.1	± 3.0
PLT	$10^3/\mu\text{L}$	70.6	± 30.0
MPV	fL	8.3	± 3.0
RDW-CV	%	15.6	± 3.0
LY%	%	61.65	± 10.00
MO%	%	3.12	± 5.00
GR%	%	35.23	± 10.00
LY	$10^3/\mu\text{L}$	1.23	± 0.50
MO	$10^3/\mu\text{L}$	0.06	± 0.20
GR	$10^3/\mu\text{L}$	0.71	± 0.50
CRP	Precision	$\leq 10.0\%$	
		(1.50-8.10 mg/L)	

Manual Registration

Lot	258
Level	Low
Control	MK-3C


<https://www.nihonkohden.com/spi/index.html>

NIHON KOHDEN CORPORATION

1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN


0694-901765K