

MEK-CAL

HEMATOLOGY CALIBRATOR

2.0mL

ASSAY SHEET

LOT	PLUS259
	2025-10-05
OPEN	7DAYS

MEK-9100,9200
2D Barcode Entry

MEK-1301,1302,1303,1305
Barcode Entry

<https://www.nihonkohden.com/spi/index.html>

NIHON KOHDEN CORPORATION

1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN


MEK-CAL



PLUS259



2025-10-05



0123



0694-901832D

		MEK -6400 -6410 -6420 -6450	MEK -6500 -6510 -6550	MEK -7300	MEK -8222	MEK -9100 -9200	MEK -1301 -1302 -1303 -1305
WBC	10 ³ /μL	9.3	9.3	9.4	9.3	9.66	9.33
RBC	10 ⁶ /μL	4.63	4.63	4.64	4.60	4.60	4.61
HGB	g/dL	13.2	13.4	13.2	13.2	13.61	13.48
HCT	%	39.7	39.7	40.5	39.4	45.4	45.6
MCV	fL	85.8	85.8	87.2	85.6	98.8	99.0
MCH	pg	28.5	28.9	28.4	28.7	29.6	29.2
MCHC	g/dL	33.2	33.8	32.6	33.5	30.0	29.6
PLT	10 ³ /μL	231	231	237	236	244.4	252.8
MPV	fL	5.6	8.3	7.1	5.6	8.3	8.4
RDW-CV *	%	16.2	16.2	16.4	15.3	16.3	16.2

MEK	
	-7300
	-8222
	-9100
	-9200
FS	172
FL	84
SD	113
FS THR.	100

* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450,8222.