

MEK-CAL

HEMATOLOGY CALIBRATOR

2.0mL

ASSAY SHEET

LOT	PLUS25X
	2025-11-05
OPEN	7DAYS

MEK-9100,9200
2D Barcode Entry

MEK-1301,1302,1303,1305
Barcode Entry

<https://www.nihonkohden.com/spi/index.html>

NIHON KOHDEN CORPORATION

1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN


MEK-CAL



PLUS25X



2025-11-05



0123



0694-901832E

		MEK -6400 -6410 -6420 -6450	MEK -6500 -6510 -6550	MEK -7300	MEK -8222	MEK -9100 -9200	MEK -1301 -1302 -1303 -1305
WBC	10 ³ /μL	8.9	8.9	8.9	8.8	9.12	8.98
RBC	10 ⁶ /μL	4.60	4.60	4.59	4.55	4.55	4.55
HGB	g/dL	13.2	13.4	13.2	13.1	13.45	13.40
HCT	%	39.7	39.7	40.2	39.1	44.7	44.9
MCV	fL	86.4	86.4	87.6	85.9	98.2	98.7
MCH	pg	28.7	29.1	28.8	28.8	29.6	29.5
MCHC	g/dL	33.2	33.8	32.8	33.5	30.1	29.8
PLT	10 ³ /μL	229	229	227	227	244.8	251.7
MPV	fL	6.3	8.6	7.4	6.3	8.8	8.8
RDW-CV *	%	16.5	16.5	16.6	15.5	16.6	16.5

MEK	
	-7300
	-8222
	-9100
	-9200
FS	167
FL	88
SD	105
FS THR.	100

* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450,8222.