

MK-3CH HEMATOLOGY CONTROL+CRP

2.5mL

ASSAY SHEET

LOT	MK3C25YH
	2026-02-05
OPEN	14DAYS

MEK-1303 Barcode Entry

If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.

1



2



3



		MEK-1303	RANGE
WBC	10 ³ /μL	20.88	± 2.20
RBC	10 ⁶ /μL	5.88	± 0.25
HGB	g/dL	19.38	± 0.70
HCT	%	62.7	± 4.0
MCV	fL	106.6	± 7.0
MCH	pg	33.0	± 2.8
MCHC	g/dL	30.9	± 3.0
PLT	10 ³ /μL	535.4	± 70.0
MPV	fL	8.3	± 3.0
RDW-CV	%	17.4	± 3.0
LY%	%	12.38	± 10.00
MO%	%	1.78	± 5.00
GR%	%	85.84	± 10.00
LY	10 ³ /μL	2.58	± 1.80
MO	10 ³ /μL	0.37	± 0.70
GR	10 ³ /μL	17.93	± 2.00
CRP	Precision	≤ 10.0% (46.80-86.80 mg/L)	

Manual Registration

Lot	25Y
Level	High
Control	MK-3C


<https://www.nihonkohden.com/spi/index.html>

NIHON KOHDEN CORPORATION

1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN

CE
0123

IVD
0694-901767L