

MEK-3DN

HEMATOLOGY CONTROL

2.0mL

ASSAY SHEET

LOT	B262N
	2026-05-05
OPEN	14DAYS

MEK-1301,1302,1305 Barcode Entry

The barcodes may be read in any order.
If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.

1



2



3



		RANGE	MEK -6400 -6410 -6420 -6450	MEK -6500 -6510 -6550	MEK -1301 -1302 -1305
WBC	10 ³ /μL	±0.8	7.6	7.6	7.61
RBC	10 ⁶ /μL	±0.20	4.80	4.80	4.75
HGB	g/dL	±0.5	13.6	13.7	13.85
HCT	%	±3.0	41.0	41.0	46.3
MCV	fL	±7.0	85.5	85.5	97.5
MCH	pg	±2.4	28.3	28.5	29.2
MCHC	g/dL	±3.0	33.2	33.4	29.9
PLT	10 ³ /μL	±50	236	236	249.9
MPV	fL	±2.0	5.7	8.2	8.2
RDW-CV *	%	±3.0	16.7	16.7	17.0
LY%	%	±10.0	27.7	27.7	31.24
MO%	%	±5.0	2.6	2.6	2.50
GR%	%	±10.0	69.7	69.7	66.26
LY	10 ³ /μL	±1.0	2.1	2.1	2.38
MO	10 ³ /μL	±0.4	0.2	0.2	0.19
GR	10 ³ /μL	±1.0	5.3	5.3	5.04

* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450.

MEK-1305
QC ESR
10±5 mm

Manual Registration

Lot	262
Level	Normal
Control	MEK-3D


<https://www.nihonkohden.com/spi/index.html>

NIHON KOHDEN CORPORATION

1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN


MEK-3DN



B262N



2026-05-05



0123



0694-901685W