

MEK-5DN
HEMATOLOGY CONTROL
for WBC 5 part differential

3.0mL

ASSAY SHEET

LOT	MEK5D262N
	2026-05-05
OPEN	14DAYS

MEK-9100,9200
2D Barcode Entry


		RANGE	MEK -7300	MEK -8222	MEK -9100 -9200
WBC	10 ³ /μL	±0.8	8.4	8.3	8.49
NE%	%	±15.0	61.5	61.5	60.83
LY%	%	±15.0	23.0	23.0	26.48
MO%	%	±5.0	8.5	8.5	6.07
EO%	%	±5.0	5.4	5.4	5.26
BA%	%	±4.0	1.6	1.6	1.36
IG% **	%	±4.0	-	-	0.48
NE	10 ³ /μL	±1.5	5.2	5.2	5.15
LY	10 ³ /μL	±1.5	1.9	1.9	2.25
MO	10 ³ /μL	±0.5	0.7	0.7	0.52
EO	10 ³ /μL	±0.5	0.5	0.4	0.45
BA	10 ³ /μL	±0.4	0.1	0.1	0.12
IG **	10 ³ /μL	±0.4	-	-	0.04
RBC	10 ⁶ /μL	±0.20	4.43	4.48	4.45
HGB	g/dL	±0.5	12.9	12.7	13.22
HCT	%	±3.0	39.4	38.7	44.9
MCV	fL	±6.0	88.9	86.4	100.8
MCH	pg	±2.8	29.1	28.3	29.7
MCHC	g/dL	±3.0	32.7	32.8	29.4
RDW-CV *	%	±4.0	15.5	14.7	15.1
PLT	10 ³ /μL	±50	284	261	257.6
MPV	fL	±2.0	5.2	3.6	6.7

Manual Registration	
Lot	262
Level	Normal
Control	MEK5D

* RDW-CV is displayed as RDW on MEK-8222.

** MEK-9200 only
Only software versions 02-01 and later can be used.



<https://www.nihonkohden.com/spi/index.html>


NIHON KOHDEN CORPORATION

1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN


MEK-5DN



262N



2026-05-05



0123



0694-901691V