

MEK-CAL

HEMATOLOGY CALIBRATOR

2.0mL

ASSAY SHEET

	PLUS262
	2026-03-05
	7DAYS

MEK-9100,9200
2D Barcode Entry


		MEK -6400 -6410 -6420 -6450	MEK -6500 -6510 -6550	MEK -7300	MEK -8222	MEK -9100 -9200	MEK -1301 -1302 -1303 -1305
WBC	10 ³ /μL	8.9	8.9	9.0	8.8	9.23	8.96
RBC	10 ⁶ /μL	4.58	4.58	4.55	4.51	4.51	4.55
HGB	g/dL	12.8	13.0	12.9	12.8	13.17	13.12
HCT	%	39.0	39.0	39.1	38.0	43.9	44.5
MCV	fL	85.2	85.2	85.9	84.3	97.3	97.8
MCH	pg	27.9	28.4	28.4	28.4	29.2	28.8
MCHC	g/dL	32.8	33.3	33.0	33.7	30.0	29.5
PLT	10 ³ /μL	245	245	239	238	252.3	264.6
MPV	fL	6.8	8.8	7.6	6.8	9.0	9.1
RDW-CV *	%	16.1	16.1	16.3	15.3	16.3	16.2

MEK	
	-7300
	-8222
	-9100
	-9200
FS	166
FL	86
SD	102
FS THR.	100

MEK-1301,1302,1303,1305
Barcode Entry

* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450,8222.


<https://www.nihonkohden.com/spi/index.html>

NIHON KOHDEN CORPORATION

1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN


MEK-CAL



PLUS262



2026-03-05



0694-901832I