

MK-3CL HEMATOLOGY CONTROL+CRP 2.5mL

ASSAY SHEET

LOT	MK3C262L
	2026-05-05
OPEN	14DAYS

MEK-1303 Barcode Entry

If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.

1



2



3



		MEK-1303	RANGE
WBC	10 ³ /μL	2.18	± 0.60
RBC	10 ⁶ /μL	2.32	± 0.15
HGB	g/dL	6.04	± 0.40
HCT	%	20.1	± 2.0
MCV	fL	86.8	± 9.0
MCH	pg	26.0	± 2.4
MCHC	g/dL	30.0	± 3.0
PLT	10 ³ /μL	74.7	± 30.0
MPV	fL	8.2	± 3.0
RDW-CV	%	16.7	± 3.0
LY%	%	61.27	± 10.00
MO%	%	3.26	± 5.00
GR%	%	35.47	± 10.00
LY	10 ³ /μL	1.34	± 0.50
MO	10 ³ /μL	0.07	± 0.20
GR	10 ³ /μL	0.77	± 0.50
CRP	Precision	≤ 10.0%	
		(2.40-13.40 mg/L)	

Manual Registration

Lot	262
Level	Low
Control	MK-3C


<https://www.nihonkohden.com/spi/index.html>

NIHON KOHDEN CORPORATION

1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN

CE
0123

IVD
0694-901765M