

**MEK-CAL** HEMATOLOGY CALIBRATOR 2.0mL

**ASSAY SHEET**

<b>LOT</b>	PLUS264
	2026-05-05
<b>OPEN</b>	7DAYS

		MEK -6400 -6410 -6420 -6450	MEK -6500 -6510 -6550	MEK -7300	MEK -8222	MEK -9100 -9200	MEK -1301 -1302 -1303 -1305
WBC	10 <sup>3</sup> /μL	8.8	8.8	8.9	8.8	9.14	8.97
RBC	10 <sup>6</sup> /μL	4.63	4.63	4.62	4.61	4.63	4.67
HGB	g/dL	13.5	13.6	13.5	13.4	13.77	13.77
HCT	%	40.7	40.7	40.9	39.8	46.0	46.7
MCV	fL	87.8	87.8	88.5	86.4	99.4	99.9
MCH	pg	29.2	29.4	29.2	29.1	29.7	29.5
MCHC	g/dL	33.2	33.4	33.0	33.7	29.9	29.5
PLT	10 <sup>3</sup> /μL	223	223	221	221	231.2	238.5
MPV	fL	6.0	8.4	7.2	6.0	8.4	8.4
RDW-CV*	%	15.8	15.8	16.0	15.0	15.8	15.7

<b>MEK</b>	
	-7300
	-8222
	-9100
	-9200
<b>FS</b>	166
<b>FL</b>	86
<b>SD</b>	103
<b>FS THR.</b>	100

MEK-9100,9200  
2D Barcode Entry



MEK-1301,1302,1303,1305  
Barcode Entry



\* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450,8222.



<https://www.nihonkohden.com/spi/index.html>

**NIHON KOHDEN CORPORATION**  
1-31-4 NISHIOCHIAI, SHINJUKU-KU  
TOKYO 161-8560 JAPAN

