


<b>MEK-3DH</b>	HEMATOLOGY CONTROL	2.0mL
----------------	--------------------	-------

**ASSAY SHEET**

LOT	B265H
	2026-08-05
OPEN	14DAYS

<b>MEK-1301,1302,1305 Barcode Entry</b>
---

The barcodes may be read in any order.  
 If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.

1



2



3



		RANGE	MEK -6400 -6410 -6420 -6450	MEK -6500 -6510 -6550	MEK -1301 -1302 -1305
<b>WBC</b>	10 <sup>3</sup> /μL	±2.2	20.2	20.2	19.89
<b>RBC</b>	10 <sup>6</sup> /μL	±0.25	5.89	5.89	5.90
<b>HGB</b>	g/dL	±0.7	19.1	19.0	19.17
<b>HCT</b>	%	±4.0	55.0	55.0	62.8
<b>MCV</b>	fL	±9.0	93.3	93.3	106.5
<b>MCH</b>	pg	±2.8	32.4	32.3	32.5
<b>MCHC</b>	g/dL	±3.0	34.7	34.5	30.5
<b>PLT</b>	10 <sup>3</sup> /μL	±70	499	499	497.1
<b>MPV</b>	fL	±3.0	5.3	7.9	8.1
<b>RDW-CV *</b>	%	±3.0	16.3	16.3	17.8
<b>LY%</b>	%	±10.0	15.8	15.8	13.30
<b>MO%</b>	%	±5.0	1.0	1.0	1.15
<b>GR%</b>	%	±10.0	83.2	83.2	85.55
<b>LY</b>	10 <sup>3</sup> /μL	±1.8	3.2	3.2	2.65
<b>MO</b>	10 <sup>3</sup> /μL	±0.7	0.2	0.2	0.23
<b>GR</b>	10 <sup>3</sup> /μL	±2.0	16.8	16.8	17.01

Manual Registration	
Lot	265
Level	High
Control	MEK-3D

\* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450.



<https://www.nihonkohden.com/spi/index.html>


**NIHON KOHDEN CORPORATION**

 1-31-4 NISHIOCHIAI, SHINJUKU-KU  
 TOKYO 161-8560 JAPAN


MEK-3DH



B265H



2026-08-05



0694-901686X