


MK-3CH HEMATOLOGY CONTROL+CRP 2.5mL

ASSAY SHEET

LOT	MK3C265H
	2026-08-05
OPEN	14DAYS

MEK-1303 Barcode Entry

If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.

1



2



3



		MEK-1303	RANGE
WBC	10 ³ /μL	20.15	± 2.20
RBC	10 ⁶ /μL	5.66	± 0.25
HGB	g/dL	18.95	± 0.70
HCT	%	60.6	± 4.0
MCV	fL	107.0	± 9.0
MCH	pg	33.5	± 2.8
MCHC	g/dL	31.3	± 3.0
PLT	10 ³ /μL	497.0	± 70.0
MPV	fL	8.0	± 3.0
RDW-CV	%	17.2	± 3.0
LY%	%	15.00	± 10.00
MO%	%	0.76	± 5.00
GR%	%	84.24	± 10.00
LY	10 ³ /μL	3.02	± 1.80
MO	10 ³ /μL	0.15	± 0.70
GR	10 ³ /μL	16.98	± 2.00
CRP	Precision	≤ 10.0%	
		(47.90-88.90 mg/L)	

Manual Registration

Lot	265
Level	High
Control	MK-3C


<https://www.nihonkohden.com/spi/index.html>

NIHON KOHDEN CORPORATION

 1-31-4 NISHIOCHIAI, SHINJUKU-KU
 TOKYO 161-8560 JAPAN


 0123


0694-901767N