


MK-3CL HEMATOLOGY CONTROL+CRP 2.5mL

ASSAY SHEET

LOT	MK3C265L
	2026-08-05
OPEN	14DAYS

MEK-1303 Barcode Entry

If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.

1



2



3



		MEK-1303	RANGE
WBC	10 ³ /μL	2.08	± 0.60
RBC	10 ⁶ /μL	2.22	± 0.15
HGB	g/dL	5.89	± 0.40
HCT	%	19.0	± 2.0
MCV	fL	85.8	± 9.0
MCH	pg	26.5	± 2.4
MCHC	g/dL	31.0	± 3.0
PLT	10 ³ /μL	67.6	± 30.0
MPV	fL	8.0	± 3.0
RDW-CV	%	16.7	± 3.0
LY%	%	66.14	± 10.00
MO%	%	3.48	± 5.00
GR%	%	30.38	± 10.00
LY	10 ³ /μL	1.38	± 0.50
MO	10 ³ /μL	0.07	± 0.20
GR	10 ³ /μL	0.63	± 0.50
CRP	Precision	≤ 10.0%	
		(2.50-14.50 mg/L)	

Manual Registration

Lot	265
Level	Low
Control	MK-3C


<https://www.nihonkohden.com/spi/index.html>

NIHON KOHDEN CORPORATION

 1-31-4 NISHIOCHIAI, SHINJUKU-KU
 TOKYO 161-8560 JAPAN


 0123


0694-901765N