


**MEK-CAL** HEMATOLOGY CALIBRATOR 2.0mL

**ASSAY SHEET**

LOT	PLUS266
	2026-07-05
OPEN	7DAYS

		MEK -6400 -6410 -6420 -6450	MEK -6500 -6510 -6550	MEK -7300	MEK -8222	MEK -9100 -9200	MEK -1301 -1302 -1303 -1305
WBC	10 <sup>3</sup> /μL	9.1	9.1	9.1	9.0	9.36	9.18
RBC	10 <sup>6</sup> /μL	4.51	4.51	4.42	4.43	4.47	4.48
HGB	g/dL	13.2	13.3	13.1	13.0	13.28	13.31
HCT	%	39.8	39.8	39.0	38.1	44.3	44.7
MCV	fL	88.2	88.2	88.3	85.9	99.2	99.7
MCH	pg	29.3	29.5	29.6	29.3	29.7	29.7
MCHC	g/dL	33.2	33.4	33.6	34.1	30.0	29.8
PLT	10 <sup>3</sup> /μL	232	232	220	219	233.6	241.3
MPV	fL	5.9	8.3	7.2	5.9	8.4	8.4
RDW-CV*	%	15.9	15.9	16.0	15.0	15.8	15.7

MEK	
	-7300
	-8222
	-9100
	-9200
FS	169
FL	86
SD	105
FS THR.	100

MEK-9100,9200  
2D Barcode Entry



MEK-1301,1302,1303,1305  
Barcode Entry



\* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450,8222.



<https://www.nihonkohden.com/spi/index.html>

**NIHON KOHDEN CORPORATION**  
1-31-4 NISHIOCHIAI, SHINJUKU-KU  
TOKYO 161-8560 JAPAN

