

MEK-CAL HEMATOLOGY CALIBRATOR 2.0mL

ASSAY SHEET

LOT	PLUS267
	2026-08-05
OPEN	7DAYS

		MEK -6400 -6410 -6420 -6450	MEK -6500 -6510 -6550	MEK -7300	MEK -8222	MEK -9100 -9200	MEK -1301 -1302 -1303 -1305
WBC	10 ³ /μL	8.8	8.8	8.8	8.7	9.02	8.93
RBC	10 ⁶ /μL	4.73	4.73	4.62	4.62	4.65	4.65
HGB	g/dL	13.6	13.8	13.5	13.4	13.76	13.77
HCT	%	40.3	40.3	40.4	39.3	45.7	45.9
MCV	fL	85.3	85.3	87.5	85.1	98.3	98.8
MCH	pg	28.8	29.2	29.2	29.0	29.6	29.6
MCHC	g/dL	33.7	34.2	33.4	34.1	30.1	30.0
PLT	10 ³ /μL	225	225	220	219	231.4	236.5
MPV	fL	5.7	8.2	7.1	5.7	8.3	8.4
RDW-CV*	%	16.2	16.2	16.2	15.1	15.9	15.9

MEK	
	-7300
	-8222
	-9100
	-9200
FS	169
FL	88
SD	106
FS THR.	100

MEK-9100,9200
2D Barcode Entry



MEK-1301,1302,1303,1305
Barcode Entry



* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450,8222.



<https://www.nihonkohden.com/spi/index.html>

NIHON KOHDEN CORPORATION
1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN

