

MEK-3DL

HEMATOLOGY CONTROL

2.0mL

ASSAY SHEET

LOT	B268L
	2026-11-05
OPEN	14DAYS

MEK-1301,1302,1305 Barcode Entry

The barcodes may be read in any order.
If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.

1



2



3



		RANGE	MEK -6400 -6410 -6420 -6450	MEK -6500 -6510 -6550	MEK -1301 -1302 -1305
WBC	10 ³ /μL	±0.6	2.0	2.0	1.96
RBC	10 ⁶ /μL	±0.15	2.39	2.39	2.32
HGB	g/dL	±0.4	6.0	6.2	6.11
HCT	%	±2.0	18.5	18.5	20.4
MCV	fL	±9.0	77.3	77.3	88.0
MCH	pg	±2.4	25.1	25.9	26.3
MCHC	g/dL	±3.0	32.4	33.5	30.0
PLT	10 ³ /μL	±30	63	63	70.1
MPV	fL	±3.0	5.8	8.3	8.2
RDW-CV *	%	±3.0	18.8	18.8	16.5
LY%	%	±10.0	56.3	56.3	58.06
MO%	%	±5.0	9.3	9.3	6.04
GR%	%	±10.0	34.4	34.4	35.90
LY	10 ³ /μL	±0.5	1.1	1.1	1.14
MO	10 ³ /μL	±0.2	0.2	0.2	0.12
GR	10 ³ /μL	±0.5	0.7	0.7	0.70

* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450.

MEK-1305
QC ESR
47±17 mm

Manual Registration

Lot	268
Level	Low
Control	MEK-3D


<https://www.nihonkohden.com/spi/index.html>

NIHON KOHDEN CORPORATION

1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN


MEK-3DL



B268L



2026-11-05



0123



0694-901684Y