

MEK-5DL
HEMATOLOGY CONTROL
for WBC 5 part differential

3.0mL

ASSAY SHEET

LOT	MEK5D268L
	2026-11-05
OPEN	14DAYS

MEK-9100,9200
2D Barcode Entry


		RANGE	MEK -7300	MEK -8222	MEK -9100 -9200
WBC	10 ³ /μL	±0.6	2.9	3.0	2.97
NE%	%	±15.0	43.7	43.7	42.41
LY%	%	±15.0	39.3	39.3	43.64
MO%	%	±5.0	10.4	10.4	7.85
EO%	%	±5.0	4.2	4.2	4.27
BA%	%	±4.0	2.4	2.4	1.83
IG% **	%	±4.0	-	-	0.62
NE	10 ³ /μL	±0.5	1.3	1.3	1.26
LY	10 ³ /μL	±0.5	1.1	1.2	1.30
MO	10 ³ /μL	±0.2	0.3	0.3	0.23
EO	10 ³ /μL	±0.2	0.1	0.1	0.13
BA	10 ³ /μL	±0.2	0.1	0.1	0.05
IG **	10 ³ /μL	±0.2	-	-	0.02
RBC	10 ⁶ /μL	±0.15	2.52	2.55	2.44
HGB	g/dL	±0.4	6.3	6.2	6.36
HCT	%	±2.0	20.3	20.0	21.8
MCV	fL	±7.0	80.5	78.4	89.2
MCH	pg	±2.4	25.0	24.3	26.1
MCHC	g/dL	±3.0	31.0	31.0	29.2
RDW-CV *	%	±4.0	17.2	15.9	16.6
PLT	10 ³ /μL	±30	45	41	38.9
MPV	fL	±3.0	5.7	5.5	7.4

Manual Registration	
Lot	268
Level	Low
Control	MEK5D

* RDW-CV is displayed as RDW on MEK-8222.

** MEK-9200 only
Only software versions 02-01 and later can be used.


<https://www.nihonkohden.com/spi/index.html>

NIHON KOHDEN CORPORATION

1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN


MEK-5DL



268L



2026-11-05



0123



0694-901690X