

MEK-5DN
HEMATOLOGY CONTROL
for WBC 5 part differential

3.0mL

ASSAY SHEET

LOT	MEK5D268N
	2026-11-05
OPEN	14DAYS

MEK-9100,9200
2D Barcode Entry


		RANGE	MEK -7300	MEK -8222	MEK -9100 -9200
WBC	10 ³ /μL	±0.8	8.0	7.9	8.16
NE%	%	±15.0	61.7	61.7	64.20
LY%	%	±15.0	23.5	23.5	24.62
MO%	%	±5.0	10.7	10.7	6.56
EO%	%	±5.0	2.9	2.9	3.47
BA%	%	±4.0	1.2	1.2	1.15
IG% **	%	±4.0	-	-	0.84
NE	10 ³ /μL	±1.5	4.9	4.9	5.24
LY	10 ³ /μL	±1.5	1.9	1.9	2.01
MO	10 ³ /μL	±0.5	0.9	0.8	0.54
EO	10 ³ /μL	±0.5	0.2	0.2	0.28
BA	10 ³ /μL	±0.4	0.1	0.1	0.09
IG **	10 ³ /μL	±0.4	-	-	0.07
RBC	10 ⁶ /μL	±0.20	4.53	4.58	4.62
HGB	g/dL	±0.5	13.5	13.4	13.74
HCT	%	±3.0	40.6	40.0	46.8
MCV	fL	±6.0	89.7	87.3	101.3
MCH	pg	±2.8	29.8	29.3	29.7
MCHC	g/dL	±3.0	33.3	33.5	29.4
RDW-CV *	%	±4.0	14.6	14.0	14.6
PLT	10 ³ /μL	±50	237	218	218.1
MPV	fL	±2.0	5.3	3.9	6.9

Manual Registration	
Lot	268
Level	Normal
Control	MEK5D

* RDW-CV is displayed as RDW on MEK-8222.

** MEK-9200 only
Only software versions 02-01 and later can be used.


<https://www.nihonkohden.com/spi/index.html>

NIHON KOHDEN CORPORATION

1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN


MEK-5DN



268N



2026-11-05



0123



0694-901691X