

MK-3CH HEMATOLOGY CONTROL+CRP

2.5mL

ASSAY SHEET

LOT	MK3C268H
	2026-11-05
OPEN	14DAYS

MEK-1303 Barcode Entry

If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.

1



2



3



		MEK-1303	RANGE
WBC	$10^3/\mu\text{L}$	20.04	± 2.20
RBC	$10^6/\mu\text{L}$	5.80	± 0.25
HGB	g/dL	19.06	± 0.70
HCT	%	61.2	± 4.0
MCV	fL	105.5	± 7.0
MCH	pg	32.9	± 2.8
MCHC	g/dL	31.1	± 3.0
PLT	$10^3/\mu\text{L}$	502.6	± 70.0
MPV	fL	8.1	± 3.0
RDW-CV	%	16.8	± 3.0
LY%	%	13.87	± 10.00
MO%	%	1.26	± 5.00
GR%	%	84.87	± 10.00
LY	$10^3/\mu\text{L}$	2.78	± 1.80
MO	$10^3/\mu\text{L}$	0.25	± 0.70
GR	$10^3/\mu\text{L}$	17.01	± 2.00
CRP	Precision	$\leq 10.0\%$ (39.60-73.40 mg/L)	

Manual Registration

Lot	268
Level	High
Control	MK-3C


<https://www.nihonkohden.com/spi/index.html>

NIHON KOHDEN CORPORATION

1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN

CE
0123

IVD

0694-9017670