

MEK-CAL

HEMATOLOGY CALIBRATOR

2.0mL

ASSAY SHEET

LOT	PLUS269
	2026-10-05
OPEN	7DAYS

MEK-9100,9200
2D Barcode Entry

MEK-1301,1302,1303,1305
Barcode Entry

<https://www.nihonkohden.com/spi>

		MEK -6400 -6410 -6420 -6450	MEK -6500 -6510 -6550	MEK -7300	MEK -8222	MEK -9100 -9200	MEK -1301 -1302 -1303 -1305
WBC	10 ³ /μL	9.0	9.0	9.1	8.9	9.30	9.10
RBC	10 ⁶ /μL	4.81	4.81	4.72	4.74	4.79	4.78
HGB	g/dL	13.9	14.0	13.8	13.7	14.04	14.04
HCT	%	41.3	41.3	40.8	39.8	46.9	46.9
MCV	fL	85.9	85.9	86.5	83.9	97.9	98.1
MCH	pg	28.9	29.1	29.2	28.9	29.3	29.4
MCHC	g/dL	33.7	33.9	33.8	34.4	29.9	29.9
PLT	10 ³ /μL	216	216	215	215	228.6	231.0
MPV	fL	5.8	8.3	7.2	5.8	8.3	8.4
RDW-CV *	%	16.1	16.1	15.9	15.0	15.9	15.6

MEK	
	-7300
	-8222
	-9100
	-9200
FS	166
FL	90
SD	105
FS THR.	100

* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450,8222.


NIHON KOHDEN CORPORATION

1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN


MEK-CAL



PLUS269



2026-10-05



0123



0694-901832P